English

Intended Use

Tuf-Link II is an addition curing vinyl silicone polymer for soft relining the hard acrylic denture base to increase patient comfort and solve difficult prosthodontic problems. Together with Dentatus Atlas Implants, Tuf-Link II is used for improved retention of overdentures. See Instructions for Atlas.

Depending on the conditions and professional preference, Tuf-Link II can be used as a chairside procedure or used indirectly in the dental laboratory. In both techniques, the most important factors are the preparation of the denture base: to reduce a uniform level of hard acrylic, and to obtain a balanced denture-tissue interface.

Instructions for Use

Prepare the denture Prepare the acrylic surface with the Tuf-Link depth control cutting tool (in Atlas system), to a depth of 1.5-2 mm. To lock in the silicone reline without adhesives, create an undercut groove encasement in the periphery of the preparation.

Gun and Cartridge Insert the cartridge in the mixing gun. Remove and discard the cartridge cap and attach the mixing tip. Bleed a small amount of material to ensure a smooth even flow.

***Note:** Store used cartridge with fixed used mixing tip. Do not re-use the cartridge cap as this may cause cross-contamination of components and thus premature curing of remaining contents. Do not store above 25°C (78°F). Do not use after expiry date. Keep away from children. For dental use only.*

Use non-latex protective gloves and a no-touch technique is recommended.

Relining the denture Apply Tuf-Link II silicone to the denture base to coat the entire cut-out surface by 2-3 mm and seat the denture in the patient’s mouth within 1 minute (at 23°C or 74°F) as material starts to set immediately.

a) Direct Chairside Application: Seat the denture and guide the patient to close the teeth into centric position. After 4-5 min. the silicone is fully set and the denture is removed from the mouth. Remove the liner from the denture, trim off excess flash material with fine scissors and replace the liner in the denture.

b) Indirect Laboratory Application: The preferred technique is to start with an impres­ sion taken with Tissue Conditioning Functional Impression Material. The impression is boxed for making a hard master model with dental stone. Mount the denture in a reline jig and prepare the denture surface as described above. Apply an alginate separator on the model before seating the Tuf-Link filled denture. Allow 12-14 minutes for the final setting of the silicone. Remove the liner from the prosthesis, cut away the excess flash material with a small sharp scissor, and snap the liner back into its indexed position.

***Note:** A spare liner can be prepared in advance for the patient emergency replacement. Subsequent relines can be made without cutting away the old material or changing the spatial relations of the prosthesis.*

Aftercare Instruct the patient to clean the prosthesis daily and avoid abrasive pastes and bleach or alcohol-containing cleaners. Patients with sufficient dexterity may be instructed to clean the liner with mild soap detergent and lukewarm tap water.

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