Dentatus Luminex® Post system

Light guides for deep cure in flared roots

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

MADE IN SWEDEN

A.

EN: LUMINEX Posts For Reinforcing Thin Walled Roots with Aesthetic Materials

1. The canal entrance is prepared in a conventional manner.
2. Insert a corresponding LUMINEX Smooth Post and mark its position with the depth indicating ring. Remove it from the canal and lubricate the post with a fine coating of vaseline.
3. Clean, acid-etch, and bond the canal in a conventional manner. Fill the canal with appropriate composite material, and reinsert the LUMINEX Post rotationally to full depth indicated by the ring stop. Polymerize the material with a light source directly and faciolingually in duration according to manufacturer's specifications. Remove the LUMINEX Post with a gentle rotational traction

Indirect Cast Post/ Core Techniques

Insert a corresponding size Smooth LUMINEX Post after the canal is reamed to obtain an impression with the post firmly indexed in the impression material.

B.

Indirect Cast Post/ Core Technique

1. After the canal is reamed to the desired size and depth, create an elongated tapered funnel at the root canal interface. Insert a corresponding size Smooth LUMINEX Post, and obtain an impression with the post firmly indexed in the impression material. Deform the Smooth LUMINEX Post with a heated spatula to a nail head shape for firm retention in the impression material
2. Lubricate post and pour a model. After the gypsum is fully set, soak the model in warm water before removing the plastic post with gentle rotational traction. Lubricate sparingly all interfaces and then insert a corresponding size Grooved LUMINEX Post for wax-up. The Grooved LUMINEX Post is incorporated in the wax-up which is burned out during the embedding procedure before casting.

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