
CePo Pilot Drills and CePo Drills – System 2

CePo Marking	Implant Ø	CePo Ø	L
=1.8	1.8 mm	1.4 mm	43.5 mm
=1.8	1.8 mm	1.4 mm	34.5 mm
=2.2	2.2 mm	1.8 mm	34.5 mm
=2.4	2.4 mm	2.0 mm	34.5 mm
=2.8	2.8 mm	2.4 mm	34.5 mm

EN English — Instructions for Use (IFU) 2	BG български — Инструкции 5
CS Český — Pokyny 8	DA Dansk — Instruktionser 11
DE Deutsch — Anleitung 14	EL Ελληνικά — Οδηγίες 17
ES Español — Instrucciones 20	ET Eesti — Kasutusjuhiseid 23
FR Français — Instructions 26	HU Magyar — Használati 29
IT Italiano — Istruzione 32	LT Lietuvis — Instrukcijos 35
LV Latviski — Norādījumus 38	NO Norsk — Instruksjoner 41
PL Polski — Instrukcji 44	PT Português — Instruções 47
RO Română — Instrucțiunile 50	SK Slovenský — Inštrukcie 53
SL Slovenščina — Navodila 56	SV Svenska — Instruktionser 59
TR Türkçe — Talimatları 62	

CePo Pilot Drills and CePo Drills: System 2 – for flapless placement

Refer to the intended implant IFU for further instructions for use: www.dentatus.com/ifu.

Material: Stainless Steel

Intended Use, Patient Group & Clinical Benefit: Creating osteotomies for dental implant installation in patients with one or more missing adult teeth.

Intended User: Licensed dentist with special training in dental implants.

Facility: Dental clinic.

Contraindications:

Absolute contraindication for dental implants:

- Patients with intravenous high dose bisphosphonate therapy.
- Patients undergoing chemo or radiotherapy.

Partial contraindication for dental implants:

- Patients who suffer from osteoporosis and medicate with oral bisphosphonates.
- Patients who have undergone radiotherapy should afterwards only be treated by specialists and only after approval of oncologist.

⚠ Warnings:

- Clean and sterilize before use.
- Reprocess if sterilization pouch is damaged before use.
- Sharp – handle with care.
- Device contaminated after use, see instructions for cleaning.

Caution:

- Drills should be used with copious amounts of biocompatible isotonic sterile solution spray to prevent overheating damage and possible bone necrosis.
- When preparing osteotomies in mandibular D1 & D2 hard bone, only sharp drills should be used.

Precautions:

- Safe handling of sharps reduces risk of infection and cross-contamination.
- Use medical gloves when handling the device.
- Extreme care should be observed to prevent accidental swallowing or aspiration. Preventive practices should always be utilized. If such an accident should occur, immediately contact a physician.
- Limitations on use: discard when dull.

Safe Disposal: Dispose drills in a sharps container to avoid accidental cuts and/or stick.

NOTE: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Cleaning, Disinfection & Sterilization:

Minimum drinking water quality should be used for cleaning and distilled or demineralized water for ultrasonic bath and final rinse. Hydrogen peroxide should not be used as a disinfectant.

Manual Cleaning:

Before first use and for reprocessing, immerse instruments in freshly prepared cleaning solution (enzymatic or low alkaline (pH≤8) detergent) as per agent manufacturer's instructions for concentration at maximum 40°C. Clean mechanically with a soft brush, working beneath the liquid level, until visibly clean. Rinse thoroughly. Clean in an ultrasonic bath for a minimum of 5 minutes at maximum 60°C, using a frequency of 35-45KHz and power min. 150W. Rinse thoroughly for min. 30 seconds.

Manual Disinfection:

Immerse into disinfection solution suitable for instruments of stainless steel. Follow manufacturer's instructions of disinfection agent for concentration and time. Disinfection validated with 70% Ethanol, 10 minutes soaking followed by air drying in safety cabinet.

Automated Cleaning and Disinfection:

Washer-disinfector (validated according to EN ISO 15883). Enzymatic or low alkaline, suited for medical devices. Non-corrosive, neutral rinsing agent, suited for medical devices. Follow agent manufacturer's instructions regarding concentration and temperature.


1. Load the instruments/disassembled parts in the washer-disinfector. Use a suitable instrument tray.
2. Run program suited for surgical instruments. Validation performed with the following parameters:
 - Pre-wash in cold water, 2x 2 min.
 - Main wash with detergent at minimum 55°C, 10 min.
 - Rinse in warm water, 2x 1 min.
 - Final rinse/disinfection in de-mineralized water at 90°C, minimum 1 min.
 - Dry at 110° C, minimum 15 min.
3. When unloading, confirm that the instruments are clean. If necessary, repeat from step 1 or use manual cleaning.

Sterilization: The instrument must be cleaned and disinfected before sterilization.

1. Place in individual sterilization pouches.
2. Check that the bag is not stretched.

3. Check correct sealing. Place pouches plastic towards plastic and paper towards paper.
4. Steam autoclave (validated according to EN 13060, EN 285, EN 17665). Run minimal cycle:
- Steam temperature/pressure: minimum 134°C (273°F) / 3.06 bar (27 psi).
 - Steam exposure time: minimum 3 min.
 - Vacuum drying: minimum 6 min.
 - Ensure the instruments are dry before storage. The instruments should be stored free from dust at 22 ± 3°C, and not exposed to direct sunlight or humidity above 70%.

	CE mark		Medical Device
	Read Instructions for use		Unique Device Identification
	Ultrasonic bath		Catalog number
	Washer-disinfector for thermal disinfection		Batch number
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified		Date of manufacture
	Right revolution (clockwise rotation) at the speed specified		Manufacturer
	Implantology		Caution

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Dentatus AB