Implants


I M P L A N T	D I A M E T E R				L E N G T H
	1.8	2.2	2.4	2.8	
R	AT-18L-R	AT-22L-R	AT-24L-R	AT-28L-R	14
M	AT-18M-R	AT-22M-R	AT-24M-R	AT-28M-R	10
O	AT-18S-R	AT-22S-R	AT-24S-R	-	7

Accessories


Tuf-Link

COMPONENTS

STEP BY STEP ILLUSTRATIONS






n-k-IFU-07_2017-04-19 © Dentatus

INSTRUCTIONS FOR USE

EN, SE, DE, IT, FR,
ES, PL, NL, LT

Dentatus Atlas® Implant System, approved for long-term use, provide stable anchor support and stabilization for problematic dentures. The Atlas Tuf-Link® lined overdentures create cushioned chewing impact and firm retention without the need of the restrictive metal coping O-ring technologies.

Dentatus denture instruments are designed for uniform depth reduction and for creating an inverted seam line encasement that holds the reline firmly in place without adhesives. The single short chairside procedure of installing the resilient liner prevents stress and tissue abrasions especially over uneven and bulbous ridges without follow-up tissue adjustment appointments.

The Atlas Dome Keeper with Tuf-Link can provide comfort and retention for full dentures and partial removables without the need of clasp or bar-clip retentive devices.

The Dentatus Narrow-Body Implant surgical restorative protocols were developed in conjunction with the Department of Implant Dentistry, New York University College of Dentistry.

Bernard Weissman

- President

CAUTIONS AND NOTES

CAUTION - Do not use an implant if the packaging is damaged or opened. The implants are intended for single use, thereby avoiding risk of infectious cross contamination if re-used.

Atlas implants, attached with the manual carrier drivers, **are delivered sterile** are designed for direct, safe, initial implant placement.

Accessories, including drills and reamers, are delivered factory clean and after removing their wrappings, they **should be cleaned and sterilized** in compliance with surgical protocols before use.

All devices with the exception of pre-sterilized implants, used in the clinical procedures should be autoclaved with steam sterilization prior to use. The following parameters have been validated in accordance with FDA recognized standard ANSI/AAMI/ISO 17665-1:2006, to provide a Sterility Assurance of 10^{-6} . Temp: 132°C, Pressure: 30psi, Time: 8 minutes. However, the health facility should validate its own autoclave steam sterilization machine in accordance with a recognized standard (in the USA: FDA approved standard).

CAUTION - Extreme care should be observed to prevent accidental swallowing or aspiration of implants or other related small accessory components used in these procedures. When possible, tools should be

attached with a fine cord to the patient's garment.

CAUTION - The CePo® Drills and reamers should be used with copious amounts of biocompatible sterile water or saline solution spray to prevent overheating damage and possible bone necrosis. When preparing osteotomies in mandibular D-1 and D-2 hard bone, only sharp drills and reamers should be used. As a safety precaution, used dull drills should be discarded.

NOTE - Special consideration should be given for the use of Atlas Implants for upper dentures. The softer maxilla bone and anatomical restrictions must be carefully considered before their use.

NOTE - These instructions inform the user about the recommended procedures when using Dentatus implant systems. They are intended for use by clinicians with at least a basic level of surgical and in-clinic implant training. It is the responsibility of the clinician to stay informed, educated and trained about general dentistry implant treatment techniques.

The printed guidelines, including Cautions and Notes, are to be regarded as additions to accepted clinical procedures and surgical protocols. Federal law restricts the sale of this device to a licensed physician or dentist.

Dentatus Narrow Body Ti alloy implants are FDA designated for long-term use and for any length of time as decided by the healthcare provider.

ATLAS CLINICAL PROCEDURES

To determine the ideal position and angle for implant placement, appropriate CT (CAT) scans, X-rays and measurements should be reviewed. Right-angle calipers can be used to measure the width of bone with special attention to the irregularities in bucco-lingual jaw concavities.

The level of bone may be estimated by using a ruler to establish the total ridge height by placing the ruler under the chin up to the crest of the ridge. An approximate volume of 6-8 mm soft tissue under the chin and over the ridge is deducted, and the remaining solid bone measurement is used as a guide to select an appropriate length implant.

A.

The locations for the implants should be indicated with a sterile marking pen on the mucosa. 4 Atlas Implants should be placed anterior to the mental foramina and separated equidistantly for a balanced, cross-arch stabilization and retention.

The Dentatus CePo Pilot 1.4 Drill is initially used to puncture the tissue (without rotation) indenting the bone for precise location of the osteotomies. The drills, operated at 1300-1500 RPM with copious amounts of biocompatible sterile water or saline solu-

tion spray, are angled for implant placement in bone. The osteotomies are prepared to the depth for the implant platform to seat in the mucosa.

The CePo Drill laser marks are designed for implant placement performed with surgical flap procedures of exposed bone. In flapless procedures, the osteotomies are made 4-5 mm deeper to compensate for the variable levels of mucosa and to prevent compression necrosis of the bone.

B.

Note: Local anesthetic infiltration is considered adequate for these procedures.

All osteotomies are initiated with the 1.4 CePo Pilot Drill. For placing implants wider than 1.8mm diameter, the openings are enlarged with the correspondingly marked drills or reamers in a straight up and down motion to the pre-drilled depth.

Note: Aligning the drill position for the osteotomies, the thumb and index finger may be placed on either side of the ridge to facilitate tactile feedback when the drill is in motion in close proximity to the mucosa.

See next page for drill systems.

Dentatus CePo Surgical Drill systems

Dentatus offers two different laser marking systems on the CePo Pilot and CePo drills. This information will guide you on how to use them.

System 2

A user-friendly marking system, which refers to the actual placement depth through the gingiva, when performing flapless placement of the implants.

The markings correspond to the 3 different lengths of Dentatus implants: short (S), medium (M) and long (L). The upper border line of each marking is used as reference to the gingival surface and is placed 11, 14 and 18 mm from tip.

The width of each marking is 2 mm.


System 1

The laser markings refer to the threaded length of the implants.

Markings are on 7, 10 and 14 mm from the tip of the implant, where the threads on the various implants end.

Marking width is approx 0.5 mm.

The implant platforms on Anew implants shall always be flush with the gingival surface.

Cat. nr.	Ø mm
MR - 1141 CePo Drill - pilot	1.8, Long
MR - 1142 CePo Drill - pilot	1.8, Short
MR - 1145 CePo Drill	2.2
MR - 1146 CePo Drill	2.4
MR - 1147 CePo Drill	2.8

Cat. nr.	Ø mm
MR - 1041 CePo Drill - pilot	1.8, Long
MR - 1042 CePo Drill - pilot	1.8, Short
MR - 1045 CePo Drill	2.2
MR - 1046 CePo Drill	2.4

C. **Implant Placement**


The sterile manual carrier driver attached to the sterilized implant is used to safely place the implant by a few threads.

The carrier driver is removed and the placement may be continued with the R/A Hpc driver (at no more than 50 RPM) or with the Ratchet, short of the platform contacting the mucosal tissue. *To obtain good initial stability, an installation torque of 40 Ncm is recommended. The R/A handpiece driver is not to be used at torque over 60 Ncm.* The placement is completed with the manual

driver that provides tactile feedback of its proper seating.

Important: make sure the implant is properly seated in the various drivers and that the driver is in a straight line with the implant; with the manual driver align the wings parallel to the slot in the implant head.

When excessive resistance is perceived, the implant is removed and the osteotomy is made slightly deeper. This will prevent high stress and possible damage to mandibular hard bone.

Assembly of Manual Adaptor and R/A Driver

1. Locate the index dot on the adapter
2. Line up the flat side of the shaft on the R/A Driver with the dot on the adapter
3. Push the R/A Driver into the opening until stop
4. Make sure the shaft end of the R/A Driver is flush with the back side of the Manual adapter

After use, remove the R/A Driver from the Adapter by pulling them apart for separate cleaning and disinfection. Clean the opening with a small interdental brush. Adapter is made in anodized aluminium and comes with an o-ring inside, to provide a better grip on the R/A Driver. When o-ring is damaged, the Adapter should be discarded.

To facilitate assembly of Adapter - R/A Driver, a minute drop of dental handpiece lubricant oil can be applied in the opening before autoclaving. Wipe off any excess oil before packaging into sterilization bag.

D.

Clinical Prosthetic Procedures

When an existing denture cannot be located spatially on the ridge aligned with the opposing teeth, the denture should be first relined with quick auto-cure resin before initiating the retrofitting procedure.

A marking pen is used to mark the outline in the denture for a seam line position that should be located 2-3 mm below the flange extensions, 1mm undercut, 4mm clearance from implants. The marking is extended up to the retromolar pads in a curved form for seamless liner adaptation.

- 1.** The half-sphere cutting instrument is used to create an inverted seam line encasement for retention of the liner without adhesives. The instrument shank is leaned toward the opposing flange, angled at 45° to create a well delineated, continuous seam line.
- 2.** The grooved instrument is used to score depth tracks in the denture without causing damage to the prepared seam line.
- 3.** The spherical cutting instrument is used to reduce the base to the depth marked levels and to create a smooth surface.
- 4.** The marking caps, placed on the Atlas Implants are used to mark their position in the denture. The caps are painted with any safe vegetable-based

colorant or marker and the denture is firmly pressed down to transfer the marked positions for the implant head.

- 5.** The spherical cutting instrument is used to indent the marked spaces allowing the denture to seat passively in place. To prevent excessive denture reduction, the instrument should be sunk down 2/3 of its spherical depth at the places of the Implants. When the marker paint no longer is transferred from the marking caps to the denture base, the reduction is sufficient.

Note: A confirming check can be performed by applying markings directly to the implant heads for a final touch up and rounding off all sharp projections to prevent incipient tears in the silicone liner.

E.

Silicone Relining Procedures

Note: The patient's mouth, tongue and denture are thoroughly rinsed to remove fine particles that may become embedded in the liner. Preparing for the quick relining procedure, the patient should be cautioned that they will need to gently press the denture in place in contact with their opposing teeth and keep firmly in place for 4-6 minutes without tongue or jaw movements.

- 1.** The Tuf-Link dispenser with an auto-mixing tip is bled before injecting the silicone into the denture base, the flow is guided into the seam line. Sufficient silicone should be used to create pressure required for a dense, nonporous liner.
- 2.** The denture is inserted from the posterior direction and guided forward with pulsating side motions. The patient is guided to contact the opposing teeth and hold still in place. The procedure must be carried out expeditiously, as the silicone begins to set within minutes.
- 3.** When polymerized, the denture is removed and the liner is gently lifted out of the denture base with a pointed instrument at one of the distal extensions.
- 4.** Check that the silicone thickness is at minimum 2 mm above the implant heads. If not, adjust the denture base accordingly.
- 5.** The excess liner is cut away along the well delineated seam line with long-nosed fine scissors.
- 6.** To reinsert the liner, it is first located at one extreme end of the denture and guided into place with light finger pressure. A smooth flat instrument may be used to guide the reline into

the seam line which will prevent fine food particle accumulation.

Atlas Tuf-Link self-aligning dentures tolerate substantial implant angulations. The denture, removed without implant stress or tissue abrasions, rarely requires follow-up tissue adjustment appointments. Depending on variable conditions, the liners should remain resilient for 12 - 18 months. Although they are quickly replaceable, a spare liner can be made for patients to keep it in a water filled container for emergency use.

POST-OP CHECKLIST

- a.** Muscle impinging areas, the frequent cause of denture dislocation in functional and parafunctional movements should be eliminated.
- b.** Occlusal contacts and tripping interferences should be corrected.
- c.** The patient should immediately report any irregularities or discomfort.

PATIENT HOME-CARE INSTRUCTIONS

Patients can immediately ingest food; however, they should observe a 48 hour soft diet regimen, without hard or chewy foods.

- a.** Patient should not remove their denture during the first two days.
- b.** Patient should not remove their denture at night for the first two weeks following Atlas placement. This will prevent any edema initiated by the surgical placement from forcing the soft tissue up, possibly interfering with the reseating of the denture on the implants. After the two week period patients may sleep with or without their Dentures.
- c.** Patient should remove denture briefly each day for oral hygiene. Instruct them to use a soft nylon bristle toothbrush with toothpaste or water to scrub the denture. Do not have patient attempt to remove the silicone liner for at least 3 days. After this initial period, patients with good dexterity may remove the liner (with a toothpick) and clean it gently in the same manner as the denture. After the initial healing period, dentures may be soaked for 30 minutes each day in a denture cleaning solution. (Note: instruct patients not to clean their dentures or silicone relines over an open sink drain). Instruct patients how to clean

the head, collar and platform of the Atlas implants.

- d.** Patient may use a pain-relieving medication to alleviate any discomfort.
- e.** Patient should be instructed to call the office if any serious pain or other symptoms arise.

RETROFITTING MAXILLA IMPLANT-SUPPORTED OVERDENTURES AND REMOVABLE PARTIALS

Retrofitting technical procedures for implant-supported dentures are substantially the same as described in the Atlas manual for mandible overdentures.

The softer, porous maxillary bone may not be able to withstand the functional forces without a splinted bar interconnecting 6-8 implants. The Dentatus Anew Implants with components design for fixed restoration can be used to create a resin-metal splint. Retrofitting the denture with the Tuf-Link silicone liner provides cushioned, comfortable, and a secure retention with a modified palatal form.

The Dentatus Anew Narrow Body Implant System, designed for reconstructive procedures are used for replacing single teeth, multi-unit restorations and full arch screw-retained fixed bridges. The minimally invasive less demanding surgical and clinical procedures are accessible and more affordable to the much larger range of patients' needs.

Careful consideration should be given to the maxilla and anatomical limitations and bone by using 6-8 implants with shallow osteotomies to half the depth made in hard bone. The self threading Anew Implants placed with the manual driver will create their own condensed thread form in the bone.

The Anew Implant illustrated technical manual for replacing missing single teeth, multi-unit and full-arch restorations are available upon request or online at www.dentatus.com.

TUF-LINK® SILICONE MATERIAL FOR USE IN THE GENERAL PRACTICE

The Dentatus Tuf-Link resilient liner that maintain stable forms over lengthy periods of time may be used to relines upper dentures, partial removables, & during sinus lift and hard and soft reconstructive procedures. The liner can also be used to line recording base forms to establish accurate occlusal, vertical and spatial relations preventing damage to master casts.

The Tuf-Link 50 ml dispenser with auto mixing tips is convenient for large volume relines. The Tuf-Link small 10 ml dispenser, is designed for small relines and to quickly relieve tissue abrasions. Two small indents made with round burs spaced 3 - 4 mm apart will securely retain the patch in place. The space may be permanently re-filled with hard resin after healing.

Other applications are in surgical procedures where tissue sites need to be safely protected over lengthy periods of time. Tuf-Link silicone material is dense and nonporous, preventing food impaction and bacterial accumulation.

ATLAS MARKETING PROGRAM

The Atlas Denture Comfort program with print and media materials can help you jump-start your implant practice to the very large number of patients that suffer with inadequate restorations. The marketing program is based around the consumer website www.denturecomfort.com featuring a Dentist locator by zip code. The program includes marketing tools such as patient brochures, posters, post-cards, press releases and much more. Call Dentatus at +1 800-323-3136 to see how you can take advantage of these opportunities, many of them free of charge (US only. For other markets, please contact Dentatus AB, www.dentatus.com).

You can also reach us through e-mail at dentatus@dentatus.com or see the full line of Dentatus distinctive products on the web at www.dentatus.com.

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