
CePo Pilot Drills and CePo Drills – System 1

| CePo Marking | Implant Ø | CePo Ø | L |
|--------------|-----------|--------|---------|
| =1.8 | 1.8 mm | 1.4 mm | 43.5 mm |
| =1.8 | 1.8 mm | 1.4 mm | 34.5 mm |
| =2.2 | 2.2 mm | 1.8 mm | 34.5 mm |
| =2.4 | 2.4 mm | 2.0 mm | 34.5 mm |
| =2.8 | 2.8 mm | 2.4 mm | 34.5 mm |

| | |
|--|---------------------------------------|
| EN English — Instructions for Use (IFU) 2 | BG български — Инструкции 5 |
| CS Český — Pokyny 8 | DA Dansk — Instruktionser 11 |
| DE Deutsch — Anleitung 14 | EL Ελληνικά — Οδηγίες 17 |
| ES Español — Instrucciones 20 | ET Eesti — Kasutusjuhiseid 23 |
| FR Français — Instructions 26 | HU Magyar — Használati 29 |
| IT Italiano — Istruzione 32 | LT Lietuvis — Instrukcijos 35 |
| LV Latviski — Norādījumus 38 | NO Norsk — Instruksjoner 41 |
| PL Polski — Instrukcji 44 | PT Português — Instruções 47 |
| RO Română — Instrucțiunile 50 | SK Slovenský — Inštrukcie 53 |
| SL Slovenščina — Navodila 56 | SV Svenska — Instruktionser 59 |
| TR Türkçe — Talimatları 62 | |

CePo Pilot Drills and CePo Drills: System 1 – for Bone Measurement

Refer to the intended implant IFU for further instructions for use: www.dentatus.com/ifu.

Material: Stainless Steel

Intended Use, Patient Group & Clinical Benefit: Creating osteotomies for dental implant installation in patients with one or more missing adult teeth.

Intended User: Licensed dentist with special training in dental implants.

Facility: Dental clinic.

Contraindications:

Absolute contraindication for dental implants:

- Patients with intravenous high dose bisphosphonate therapy.
- Patients undergoing chemo or radiotherapy.

Partial contraindication for dental implants:

- Patients who suffer from osteoporosis and medicate with oral bisphosphonates.
- Patients who have undergone radiotherapy should afterwards only be treated by specialists and only after approval of oncologist.

⚠ Warnings:

- Clean and sterilize before use.
- Reprocess if sterilization pouch is damaged before use.
- Sharp – handle with care.
- Device contaminated after use, see instructions for cleaning.

Caution:

- Drills should be used with copious amounts of biocompatible isotonic sterile solution spray to prevent overheating damage and possible bone necrosis.
- When preparing osteotomies in mandibular D1 & D2 hard bone, only sharp drills should be used.

Precautions:

- Safe handling of sharps reduces risk of infection and cross-contamination.
- Use medical gloves when handling the device.
- Extreme care should be observed to prevent accidental swallowing or aspiration. Preventive practices should always be utilized. If such an accident should occur, immediately contact a physician.
- Limitations on use: discard when dull.

Safe Disposal: Dispose drills in a sharps container to avoid accidental cuts and/or stick.

NOTE: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in

which the user and/or patient is established.

Cleaning, Disinfection & Sterilization:

Minimum drinking water quality should be used for cleaning and distilled or demineralized water for ultrasonic bath and final rinse. Hydrogen peroxide should not be used as a disinfectant.

Manual Cleaning:

Before first use and for reprocessing, immerse instruments in freshly prepared cleaning solution (enzymatic or low alkaline (pH≤8) detergent) as per agent manufacturer's instructions for concentration at maximum 40°C. Clean mechanically with a soft brush, working beneath the liquid level, until visibly clean. Rinse thoroughly. Clean in an ultrasonic bath for a minimum of 5 minutes at maximum 60°C, using a frequency of 35-45KHz and power min. 150W. Rinse thoroughly for min. 30 seconds.

Manual Disinfection:

Immerse into disinfection solution suitable for instruments of stainless steel. Follow manufacturer's instructions of disinfection agent for concentration and time. Disinfection validated with 70% Ethanol, 10 minutes soaking followed by air drying in safety cabinet.

Automated Cleaning and Disinfection:

Washer-disinfector (validated according to EN ISO 15883). Enzymatic or low alkaline, suited for medical devices. Non-corrosive, neutral rinsing agent, suited for medical devices. Follow agent manufacturer's instructions regarding concentration and temperature.

1. Load the instruments/disassembled parts in the washer-disinfector. Use a suitable instrument tray.
2. Run program suited for surgical instruments. Validation performed with the following parameters:
 - Pre-wash in cold water, 2x 2 min.
 - Main wash with detergent at minimum 55°C, 10 min.
 - Rinse in warm water, 2x 1 min.
 - Final rinse/disinfection in de-mineralized water at 90°C, minimum 1 min.
 - Dry at 110° C, minimum 15 min.
3. When unloading, confirm that the instruments are clean. If necessary, repeat from step 1 or use manual cleaning.


Sterilization: The instrument must be cleaned and disinfected before sterilization.

1. Place in individual sterilization pouches.
2. Check that the bag is not stretched.
3. Check correct sealing. Place pouches plastic towards plastic and paper towards

paper.

4. Steam autoclave (validated according to EN 13060, EN 285, EN 17665). Run minimal cycle:

- Steam temperature/pressure: minimum 134°C (273°F) / 3.06 bar (27 psi).
- Steam exposure time: minimum 3 min.
- Vacuum drying: minimum 6 min.
- Ensure the instruments are dry before storage. The instruments should be stored free from dust at 22 ± 3°C, and not exposed to direct sunlight or humidity above 70%.

| | | | |
|---|---|---|------------------------------|
|  | CE mark |  | Medical Device |
|  | Read Instructions for use |  | Unique Device Identification |
|  | Ultrasonic bath |  | Catalog number |
|  | Washer-disinfector for thermal disinfection |  | Batch number |
|  | Sterilizable in a steam sterilizer (autoclave) at the temperature specified |  | Date of manufacture |
|  | Right revolution (clockwise rotation) at the speed specified |  | Manufacturer |
|  | Implantology |  | Caution |

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Dentatus AB