
EN - Instruction for use  
DE - Gebrauchsanweisung  
IT - Istruzioni per l'uso  
FR - Mode d'emploi  
ES - Instrucciones de uso  
SE - Bruksanvisning  
NL - Instructies  
LT - Instrukcijos  
EL - Οδηγίες χρήσης

# ANEW<sup>®</sup> COMPONENTS


**Catalog no  
REF / Item**

**Product Description**

SR-49	2x ANEW Titanium Indexed Platform (for Large Spans)	
SR-96-2	2x ANEW Resin Indexed Platform (for Large Spans)	
SR-111-2	2x ANEW Titanium Indexed Platform (for Individual Crowns)	
SR-97-2	2x ANEW Resin Indexed Platform (for Individual Crowns)	
SR-189-2	2x ANEW Resin Castable Platform	
SR-58-2	2x ANEW Brass Cap Plug	
SR-68-2	2x ANEW Resin Impression Coping	
SR-123-2	2x ANEW Metal Impression Coping	
SR-62-2	2x ANEW Resin Screwcap - White	
SR-75-2	2x ANEW Resin Screwcap - Black	
SR-45-2	2x ANEW Titanium Screwcap	
SR-56	12x ANEW Gingival Spacer	
SR-83-2	2x ANEW Fixed Lab Analog	
SR-122-2	2x ANEW Removable Dowel Analog	

## Anew<sup>®</sup> Implants

		Diam (mm)	Tot. length (mm)	Thread length (mm)
ST-18S-R	ANEW Implant Pack	1.8	14	7
ST-18M-R	ANEW Implant Pack	"	17	10
ST-18L-R	ANEW Implant Pack	"	21	14
ST-22S-R	ANEW Implant Pack	2.2	14	7
ST-22M-R	ANEW Implant Pack	"	17	10
ST-22L-R	ANEW Implant Pack	"	21	14
ST-24S-R	ANEW Implant Pack	2.4	14	7
ST-24M-R	ANEW Implant Pack	"	17	10
ST-24L-R	ANEW Implant Pack	"	21	14
ST-28M-R	ANEW Implant Pack	2.8	17	10
ST-28L-R	ANEW Implant Pack	"	21	14

## Table of Contents

ENGLISH	Instruction for use	5
DEUTSCH	Gebrauchsanweisung	19
ITALIANO	Istruzioni per l'uso	33
FRANÇAIS	Mode d'emploi	47
ESPAÑOL	Instrucciones de uso	61
SVENSKA	Bruksanvisning	75
NEDERLANDS	Instructions	89
<b>LIETUVIŲ K.</b>	<b>Instrukcijos</b>	<b>103</b>
<b>ΕΛΛΗΝΙΚΑ</b>	<b>Οδηγίες χρήσης</b>	<b>117</b>

ENGLISH

DEUTSCH

ITALIANO

FRANÇAIS

ESPAÑOL

SVENSKA

NEDERLANDS

LIETUVIŲ K.

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Case Studies / Fallstudien / Caso di studio /  
Études de cas / Caso práctico / Fallstudier / Geval  
studies / **Atvejų turimai** / **Οι περιπτώσιολογικές μελέτες**

131

## Dentatus Implant Technologies – ANEW Implant System

Dear Doctor,

The ANEW Narrow Body Implant system is intended for fixed and removable dental prosthetics, and for laboratory constructed restorations. The restorations, attached to the implants with non-hygroscopic resin screwcaps without cements, can be dismantled and re-assembled without force or damage to the implants or restorations.

ANEW Narrow Body Implants are intended for patients with narrow bone, constricted inter-root and slender crown spaces without orthodontic and lengthy bone grafting interventions. Interim restorations are completed in a productive chairside procedure, replacing the patient's missing teeth.

ANEW uniform interchangeable platforms provide minimally invasive restorative options without delays and with minimal patient discomfort, as stored temporaries can be reattached in lengthier technical procedures.

The Dentatus Implant restorative protocol was developed in conjunction with the Department of Implant Dentistry, New York University College of Dentistry.

Dentatus Ti Alloy Narrow Body Implants are designated for long-term use or any length of time as decided by the health care provider.

*Bernard Weissman*  
- President

## CAUTIONS AND NOTES

**CAUTION** - Do not use an implant if the packaging is damaged or opened. The implants are intended for single use, thereby avoiding risk of infectious cross contamination if re-used.

A new implants, attached with the manual carrier drivers, are delivered sterile and are redesigned for direct, safe, initial implant placement.

Accessories, including drills, are delivered factory clean and after removing their wrappings, they should be cleaned and sterilized in compliance with surgical protocols before use.

All devices with the exception of pre-sterilized implants, used in the clinical procedures should be autoclaved with steam sterilization prior to use. The health facility should validate its own autoclave steam sterilization machine in accordance with a recognized standard (in the USA: FDA approved standard).

The following parameters have been validated in accordance with recognized standard ANSI/AAMI/ISO 17665-1:2006, to provide a Sterility Assurance of  $10^{-6}$ . Temp: 132°C, Pressure: 30psi, Time: 8 minutes.

**CAUTION** - Extreme care should be observed to prevent accidental swallowing or aspiration of implants or other related small accessory components used in these procedures. When possible, tools should be

attached with a fine cord to the patient's garment.

**CAUTION** - The CePo® Drills should be used with copious amounts of biocompatible isotonic sterile solution spray to prevent overheating damage and possible bone necrosis. When preparing osteotomies in mandibular D-1 and D-2 hard bone, only sharp drills should be used. As a safety precaution, used dull drills should be discarded.

**NOTE** - These instructions inform the user about the recommended procedures when using Dentatus implant systems. They are intended for use by clinicians with at least a basic level of surgical and in-clinic implant training. It is the responsibility of the clinician to stay informed, educated and trained about general dentistry implant treatment techniques.

The printed guidelines, including Cautions and Notes, are to be regarded as additions to accepted clinical procedures and surgical protocols. In USA Federal law restricts the sale of this device to a licensed physician or dentist.

Dentatus Narrow Body Ti alloy implants are FDA designated for long-term use or for any length of time as decided by the healthcare provider.

### Pre-Operative Planning Procedures

The patient's health history and treatment plan should be formalized to obtain the patient's informed consent before implant installation and before commencing the restorative procedures.

### Implant Site Evaluation

Radiographs are used to evaluate bone levels and height of the receptor site. Ridge anatomy and bone density can be determined with cone beam computerized tomography (CBCT) technologies.

### Flapless Protocol

The Anew Narrow Body Implants are frequently installed through the tissue without raising a flap. Single implant placement, for the support of the crown, can be aligned without stents in a confirmative manner with adjacent teeth.


When a surgical flap is raised, the procedure must be performed under strict surgical ground rules. The surgically exposed bone should be indented with a small round instrument marking the implant position. After implant installation, a primary non-tension closure is made with appropriate sutures. The suture site may be covered with a thin sterile adhesive tape to avoid suture entrapment in the resin used in the restorative procedures.

For multi-unit and full arch implants supported restorations a surgical guide is essential to establish proper occlusal alignment with esthetic crown forms.

## Protocol for Flapless Site Preparation


1. The site is swabbed with an iodine disinfectant solution or thoroughly rinsed with chlorhexidine.

2. Before installing implants in the esthetic zone the gingival tissue is prepared to create the minimum required 5mm vertical height for Screwcap assembly and for esthetic emergence profiles for crowns.

## Establishing Implant Position Angulation and Trajectory

3. The CePo 1.8mm Pilot Drill operated at 1,300-1,500RPM, with copious irrigation with biocompatible sterile water or saline solution is used to initiate the osteotomies.

**NOTE:** To establish the correct implant angulation in limiting spaces, the channel should be prepared only to half the intended depth and use x-ray to confirm its' projected trajectory. Angle corrections can be made by firmly guiding the drill to the ideal position and intended depth without materially enlarging the channels.

## Dentatus CePo Surgical Drill systems


Dentatus offers two different laser marking systems on the CePo Pilot and CePo drills. This information will guide you on how to use them.

### System 2

User-friendly marking system, which refers to the actual placement depth through the gingiva, when performing flapless placement of the implants.

The markings correspond to the 3 different lengths of Dentatus implants: short (S), medium (M) and long (L). The upper borderline of each marking is used as reference to the gingival surface and is placed 11, 14 and 18 mm from tip.

The width of each marking is 2 mm.

### System 1

The laser markings refer to the threaded length of the implants.

Markings are on 7, 10 and 14 mm from the tip of the implant, where the threads on the various implants end.

Marking width is approx 0.5 mm.

The implant platform on Anew implants shall always be flush with the gingival surface.

Cat. nr.	Ø mm
MR - 1141 CePo Drill - pilot	1.8, Long
MR - 1142 CePo Drill - pilot	1.8, Short
MR - 1145 CePo Drill	2
MR - 1146 CePo Drill	2.4
MR - 1147 CePo Drill	2

Cat. nr.	Ø mm
MR - 1041 CePo Drill - pilot	1.8, Long
MR - 1042 CePo Drill - pilot	1.8, Short
MR - 1045 CePo Drill	2.2
MR - 1046 CePo Drill	2

## Establishing Osteotomy Depth for Implants

4. The essential consideration is for the implant platform to be seated in the variable soft tissue levels. Consequently a 1-1.5mm deeper channel is made for implant full insertion.

## Enlarging the Osteotomy Channel


5. The CePo Drills correspondingly marked for 2.2 or 2.4mm Implants are used to enlarge the channels with a straight up and down pass to the pilot established depths. [\(read more about the drill systems on the previous page\)](#)

NOTE: The range of Implants and CePo Drill sizes should be available at chairside, as the inevitable variables of bone density and volume become apparent during surgery.

Having the full range of implants and drills available may prevent a compromise or the need to postpone the procedure.


## Installation of Anew Implants

1. The Anew sterilized implants with Gingival Spacers are pre-attached to the blue Carrier Driver, which is used to partially install the implant. The Carrier Driver will automatically disengage by the resistance of hard bone.

2. To obtain good initial stability, an installation torque of 40 Ncm is recommended. The R/A Handpiece Implant Driver, at a speed-controlled motor up to 50 RPM, is used with moderate downward pressure to drive the implants short of the platform touching the gingival tissue. The R/A Handpiece Driver is not to be used at torque over 60 Ncm.

3. The installation is completed with the Manual Adapter for R/A Driver, to make the upper platform level to be flush within the gingival tissue. To facilitate the crown construction, the squared sides of the platform of the implants should follow the curve of the jaw.

#### Assembly of Manual Adaptor and R/A Driver

1. Locate the index dot on the adapter
2. Line up the flat side of the shaft on the R/A Driver with the dot on the adapter
3. Push the R/A Driver into the opening until stop
4. Make sure the shaft end of the R/A Driver is flush with the back side of the Manual adapter


After use, remove the R/A Driver from the Adapter by pulling them apart for separate cleaning and disinfection. Clean the opening with a small interdental brush. Adapter is made in anodized aluminium and comes with an o-ring inside, to provide a better grip on the R/A Driver. When o-ring is damaged, the Adapter should be discarded.

To facilitate assembly of Adapter-R/A Driver, a minuted drop of dental handpiece lubricant oil can be applied in the opening before autoclaving. Wipe off any excess oil before packaging into sterilization bag.

4. The Gingival Spacer mounted on the implant is left in place during restorative procedures to prevent hard resin entrapment of the crown. The spacer is removed after the installation of the completed temporary crown.

## Chairside Construction of Temporary / Interim Crowns

1. With the Gingival Spacer tightly in place, the Resin Indexed Platform for temporary crowns is seated on the implant square platform and firmly attached with the black Screwcap Driver.


2. The Cap Plug is inserted into the Screw-cap to prevent flowing resin from blocking its access. Concave gingival interferences in adjacent teeth should be blocked out for unimpeded crown insertion.
3. A trephine drill is used for making lingual or occlusal opening in the temporary crown/form to place it over the protruding screw cap.

4. Malleable stage auto-cure resin is firmly adapted around the Resin Platform and the Screwcap assembly. The resin filled temporary crown is seated, and excess resin is removed before polymerization. Slight wrap-around contacts of the provisional/temporary crown to adjacent teeth will stabilize the implant during osseointegration and are removed at a later time.


5. The Cap Plug is removed and the Screwcap is disassembled with a Screwcap Driver.

6. Restorations can be attached to any of the labanalog with the black Screwcap. The Dowel analog can also be used as a convenient handle for touch-ups of delicate restorations.

7. Before attaching the temporary crown to the implant, the Gingival Spacer is removed from the implant and the site is cleaned and rinsed to remove impacted particles.


8. The restoration is securely assembled chairside with the white Screwcap, which is reduced slightly below the crown level. The Screwcap opening is blocked with a small dab of a cotton pellet and sealed with off color resin for easy re-access.

## Impression and technical procedures

NOTE: The resin seal in the temp-crown is removed for inserting the Chisel-End Screwcap Driver into the Screwcap to dismantle the restoration.

a) The roundside of the Impression Coping should be always oriented facially to prevent errors of assembly. The coping is firmly attached to the implant with the black Screwcap and is sealed with the Cap Plug to prevent entry of the impression material.



b) Impression Copings impinging on tissue or in close contact with adjacent teeth can be adjusted. With multiple implants in the impression, the adjusted copings must be indelibly marked for correct replacement into the impression.

c) Elastomeric impression material in a closed tray is used for full arch impressions.

d) The polymerized impression is removed with the attached coping left in place on implant.

e) The Cap Plug is removed and the Impression Coping is disassembled from the implant. The Coping, re-assembled manually to the Analog is reinserted without the Cap Plug, into the impression's indexed position.

NOTE: The impression sent to the lab should include instructions, bite registration for articulating the models, and appropriate new components, for constructing ceramo-metal or resin restorations.

Analog

## Laboratory Technical Procedures

After the model is made, the Analogs are cleaned before attaching the indexing Resin Castable Platform to the Analog with the black Screwcap.

NOTE: Made of non-hygroscopic resin, the Screwcaps require no lubricants to repel restorative material; however, when wax-up is delicate/thin the Screwcaps may be lightly coated with a suitable insulation liquid.

The wax pattern surrounding the platform and screwcaps can be cast in dental metal of choice with conventional lab procedures.

The cast metal frame with space is prepared for aesthetic ceramic or resin crown build-up.

NOTE: The castings are checked for fit on the platform and for unimpeded Screwcap insertion. Corrections made on the model or on the restoration should be clearly marked for the dentist's attention.

Chairside Assembly of Restorations  
Adjustments and reduction of hard metal should be made extraorally to prevent adverse heat effect on the implant bone interface.

The white Screwcap used for attaching restoration is reduced slightly below crown level and blocked with a small dab of a cotton pellet. It is sealed with off-color resin for re-access.

## Post-Operative Implant Care

Implants typically need a period of four to six months to integrate, depending upon the patient's healing properties and the type of bone they have. The patient must give every consideration to keep the surgical site clean and free of food particles. The following instructions are for your patient's care after surgery.

### 1. No Smoking, No Straws

Any type of suction can disturb the blood clot and loosen sutures. Smoking is to be avoided since it increases the heat in the surgical area and significantly lowers the body's ability to heal the site. Straws cannot be used for three days in most cases.

### 2. Gentle Rinsing

Frequent gentle rinsing with chlorhexidine solution for dental/oral use is recommended during a limited time. Read instruction from manufacturer of solution.

The use of commercial alcohol containing mouthwashes during the healing period is not recommended.

### 3. Ultra Soft Brush

Begin using an ultra-suave brush 24 hours after surgery to cleanse and stimulate the tissue. Gentle brushing will not harm the sutures. Brush the area at least 3 times a day. Rinse with dental chlorhexidine solution.

#### 4. Discomfort

It is very normal for patients to have some discomfort after dental procedures. If any medications are prescribed, the dosage instructions and warnings should be followed. Contact the care provider if any pain or discomfort persists.

#### 5. Swelling

Some swelling is to be expected and is not unusual. In most instances swelling can be prevented and controlled. Peak swelling is usually between 36 and 48 hours post operatively. Be advised that if there's still swelling or pain after three days, contact care provider immediately.

#### 6. Diet

Following surgery it is best to restrict the diet to fluids and soft foods for the first 24 hours. A normal diet may be resumed after this. If there is difficulty chewing, either blenderize the food or take a diet supplement. Hot foods should be avoided for the first 48 hours. Only cold/room temperature foods should be eaten.

#### 7. Implant Follow-up

Patients should see care provider for periodic check-ups until the implants become stable and integrated into the bone.

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



Anew Implant-supported Crowns Completed in a Single Chairside Procedure  
Courtesy of Dr. Sang-Choon Cho




Single tooth replaced with Anew Implant-supported ceramo-metal restoration  
Courtesy of Dr. Paul Petrungaro

Lower Anterior Crown Replaced with Anew Implant-supported  
Temporary & a Ceramo-Metal Crown Restoration  
Courtesy of Dr. Tommaso Ravasini




Lower Anterior Bridge supported by  
2 Anew Implants Constructed in a Chairside Procedure

## A 2-Unit Anew Implant Bridge Replacing Instantly a Partial Denture and Stabilizing a Bone Graft

Courtesy of Dr. Sang-Choon Cho, NYU Implant Department

### CePo® Drills & Sizing Reamers

MR-1041	2x 1.4 CePo Pilot Drills LONG (= 1.8 Impl. Diam.)
MR-1042	2x 1.4 CePo Pilot Drills SHORT (= 1.8 Impl. Diam.)
MR-1045	2x CePo Drills (= 2.2 Impl. Diam.)
MR-1046	2x CePo Drills (= 2.4 Impl. Diam.)
MR-1141	2x 1.4 CePo Pilot Drills LONG (= 1.8 Impl. Diam.) System 2
MR-1142	2x 1.4 CePo Pilot Drills SHORT (= 1.8 Impl. Diam.) System 2
MR-1145	2x CePo Drills (= 2.2 Impl. Diam.) System 2
MR-1146	2x CePo Drills (= 2.4 Impl. Diam.) System 2
MR-1147	2x CePo Drills (= 2.8 Impl. Diam.) System 2


# ANEW<sup>®</sup> COMPONENTS

**Catalog no  
REF / Item**

**Product Description**


SR-61-2	2x ANEW Titanium Connective Bar
SR-3-2	2x ANEW Short Screwcap Driver
SR-170-1	1x ANEW Long Screwcap Driver
MR-120	1x Manual Adaptor for R/A Driver
SR-57-1	1x ANEW R/A Hpc. Implant Driver
SR-30-1	1x ANEW Manual Implant Driver
SR-59-1	1x ANEW Trephine Hollow Drill
SR-5	2x ANEW Screwcap Space Refining Reamer
SR-46-1	1x ANEW Chisel-End Screwcap Driver
MPD	1x Dentatus Manual Power Driver - Handle
MPD-ST-1	1x Anew Driving Bit for MPD handle
SK-900-120	1 x Holder, Manual Adaptor in Atlas Sterilizing Box
SK-900	1x ANEW Sterilizable Instrumentation Kit Box (no instruments included)
MR-900-25	25x Moisture Absorbing Tray Papers, for SK-900

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**Dentatus AB**