
Luscent Anchors®

EN	English — Instructions	1	BG	български — Инструкции	6
CS	Český — Pokyny	11	DA	Dansk — Instruktioner	16
DE	Deutsch — Anleitung	21	EL	Ελληνικά — Οδηγίες	26
ES	Español — Instrucciones	31	ET	Eesti — Kasutusjuhiseid	36
FI	Finnish — Käyttöohjeet	41	FR	Français — Instructions	46
HR	Hrvatski — Upute za upotrebu	51	HU	Magyar — Használati	56
IT	Italiano — Istruzione	61	LT	Lietuvis — Instrukcijos	66
LV	Latviski — Norādījums	71	NO	Norsk — Instruksjoner	76
PL	Polski — Instrukcji	81	PT	Português — Instruções	86
RO	Română — Instrucțiunile	91	SK	Slovenský — Inštrukcie	96
SL	Slovenščina — Navodila	101	SV	Svenska — Instruktioner	106
TR	Türkçe — Talimatları	111			

Luscent Anchors® INSTRUCTIONS FOR USE ENGLISH

This IFU is available for download on our website: www.dentatus.com/ifu.

These instructions inform the user about the recommended procedures when using **Luscent Anchors**. They are intended for use by clinicians with a basic level of restorative dentistry and endodontic training. It is the responsibility of the clinician to stay informed, educated and trained. The printed guidelines, including Precautions and Notes, are to be regarded as additions to accepted clinical procedures and protocols.

Luscent Anchors Kits consist of reamers and specially fitted posts in different sizes. The posts are placed temporarily or permanently in the prepared root canal. The anatomically shaped reamers and posts correspond to the natural morphology of the root canal.

Dentatus Reamers are designed to be used in standard contra angle handpieces, at low speed, not exceeding 10,000 rpm. The shape of the reamers corresponds to the shape of the Luscent Anchors (see Chart 1). The Reamers are available in various lengths, and in six different diameters, #1-6. Actual diameter is shown with corresponding number of grooves on the shank (see Chart 2). The reamer system is size matched to ensure passive seating, when utilizing a size 3 post with a size 3 reamer there is no engagement of the canal walls. This prevents tension build-up and risk of root fractures.

Material:	Post: Transparent glass fiber/epoxy solid rod. <u>Accessories</u> Reamer: Stainless Steel Depth Marker: PVC-compound.	Intended Purpose:	Prefabricated fiberglass anchorage posts for retention of core materials in endodontically treated teeth.
Intended User:	Licensed Dentist.	Patient population group:	Patients with permanent teeth that are root filled and have extensive coronal damage.
Expected Clinical Benefit:	Luscent Anchors have retention properties with all standard restorative resin materials and suit a wide array of indications and requirements. They provide more esthetic results as compared to metal posts, with reliable and proven results.	Performance Characteristics:	Luscent Anchors and reamers are a state-of-the-art post system for retention of core materials in endodontically treated teeth with improved mechanical adhesion via a roughened surface. Luscent Anchors transmit light, allowing for complete polymerization of composite inside the canals and reflect natural hues for improved esthetics.
Contra-indications:	Patients suffering from bruxism or suspected bruxism, those with deep overbites, and those with insufficient crown to root ratio. Patients with known allergy to materials: reamers contain nickel.	Safe Disposal:	Posts and depth stops shall be disinfected, then the devices can be disposed of in normal waste in the clinic, according to local regulations. Dull reamers shall be disinfected, then the device can be disposed of in normal metal waste in the clinic, according to local regulations.

PRECAUTIONS:

- Delivered nonsterile and should be disinfected prior to use.
- Posts are intended for single use to avoid risk of infectious cross contamination.
- Reamers are delivered factory clean and should be cleansed and sterilized before use according to provided disinfection and sterilization instructions for an aseptic procedure.
- All other instrumentation used in the clinical procedure should be autoclaved with steam sterilization prior to use. The facility should validate its own autoclave steam sterilization machine in accordance with a recognized standard.

CAUTION:

Extreme care should be observed to prevent accidental swallowing or aspiration of endodontic posts or other related components used in this procedure. Preventive practices (rubber dam, floss ties or throat pack) should always be utilized. If such an accident should occur, immediately contact a physician.

- Dual cure adhesives and restorative materials may safely be combined with Luscent Anchors. Self-curing resins that may prematurely polymerize in the canal should be avoided.
- Damaged posts should be discarded.
- Posts are designed for passive cementation in root canals to prevent risk of root fracture.

PRE-USE INSTRUCTIONS FOR POSTS:

Posts are single use items; disinfect them with alcohol and air dry before use.

REAMERS REPROCESSING INSTRUCTIONS:

Reprocessing instructions in accordance with ISO 17664-1.

Caution:	Only use chemicals suitable for stainless steel. Reamers are delivered factory clean and after removing their wrappings they shall be cleansed, disinfected and sterilized before first use and between uses according to the instructions below.
Limitations on Reprocessing:	Reamers are intended for re-use and are delivered factory clean. The devices should be cleansed and sterilized before use according to provided disinfection and sterilization instructions for an aseptic procedure. Discard reamers when dull or damaged.
Initial Treatment at the Point of Use:	PP Instruction: Wipe off the devices after use to prevent soil and debris to dry onto the instrument. Perform cleaning as soon as possible after use. Do not exceed 2 hours.

CLEANING: MANUAL & ULTRASONIC BATH

Equipment:	Soft bristle brushes of various sizes and ultrasonic bath.
Detergent:	Follow the agent manufacturers recommendations for concentration and temperature. Manual: Enzymatic or low-alkaline (pH ≤8) detergent suitable for manual cleaning. Ultrasonic bath: Enzymatic or low-alkaline detergent with minimal foaming characteristics.

How to Use:

After endodontic therapy, root-filling material is removed to the predetermined depth with a Pathfinder, Gates-Glidden drill, peeso reamer and/or hot instrument (Fig.1). A minimum 4 mm of the root-filling material should remain apically. Radiographic verification is recommended.

The preparation should include at least 1.5 mm ferrule of sound tooth structure around the circumference of the preparation. The preparation is commenced by using the Conical [RUH series] Reamers with a slow speed contra-angle in sequential order until the desired preparation has been achieved (Fig.2). Preparation depth can be controlled by using the Depth Stop [MRA] on the reamers. In larger canals utilizing post size 4 and above, the Probos II [RUD-4] can be used to create retention grooves in the canal for additional mechanical retention against debonding.

Choosing Post Dimension:

The appropriate size post is a combination of the diameter of the reamer, the length of the canal, and the height of the crown. The sleeve on the post is color coded to correspond with the post's diameter. The grooves on the shaft of the reamer correspond with the reamer's diameter (see Chart 2).

The size number of the post coincides with the last used reamer and the length should be as long as possible without the post head interfering with shape, function, and esthetic properties of the finished restoration.

The fit in the root canal is confirmed without rotating the post and the colored sleeve is moved to the intended post height. The excess is cut with a water-cooled diamond instrument and the post are disinfected with alcohol & air-dried (Fig.4). The prepared root canal shall be acid etched then thoroughly cleansed and dried with paper points prior to placement of the post (Fig.5).

Core Build-up:

Use a flowable dental composite of choice according to the manufacturer's instructions. The bonding agent is applied to the post and in the opening of the prepared root canal to coat the wall (Fig.6). Light-cure composite is introduced into the canal (Fig.7). Dual-cure resin is recommended in canals deeper than 10mm. Use an instrument of choice (i.e.: locking tweezer) to slowly insert the post to full depth allowing excess composite to vent (Fig.8). The light cure is placed directly over the Luscent Anchor supporting it in place (Fig.9). Light cure for 1-2 minutes depending on volume and depth. The core can now be prepared for the final aesthetic restoration (Fig.10).

Post-Op:

Provide the patients with post operative instructions on hygiene and maintenance. In the event of any malfunction, patients should contact their dental provider.

Dentatus conforms to the vigilance system according to EU requirements. In the event of serious incident in relation to the device, events shall be reported to Dentatus and the competent authority of the Member State in which the provider and/or patient is established without delay.


POSTS & CORRESPONDING REAMERS

Surfer	Plain Luscent	Helix Classic Reamers	Standard Classic Reamers
		Short 28 mm Long 33 mm X Long 37 mm	Short 28 mm Long 33 mm
#1	-	RSA-1/3 RSB-1/3 RSC-1/2	RUA-1/3 RUB-1/3
#2	-	RSA-2/3 RSB-2/3 RSC-2/2	RUA-2/3 RUB-2/3
#3	White (XS)	RSA-3/3 RSB-3/3 RSC-3/2	RUA-3/3 RUB-3/3
#4	Yellow (S)	RSA-4/3 RSB-4/3 RSC-4/2	RUA-4/3 RUB-4/3
#5	Red (M)	RSA-5/3 RSB-5/3 RSC-5/2	RUA-5/3 RUB-5/3
#6	Blue (L)	RSA-6/3 RSB-6/3	RUA-6/3 RUB-6/3
	Luscent	Conical Reamers	Conical Reamers
	Purple (XXXS)		RUG-1/2
	Green (XXS)		RUG-2/2
	White (XS)	RSG-3/3 RSH-3/3 RSI-3/2	RUG-3/2
	Yellow (S)	RSG-4/3 RSH-4/3 RSI-4/2	RUG-4/2 RUH-4/2
	Red (M)	RSG-5/3 RSH-5/3 RSI-5/2	RUG-5/2 RUH-5/2
	Blue (L)	RSG-6/3 RSH-6/3	RUG-6/2 RUH-6/2

Chart 1


#1	#2	#3	#4	#5	#6	#3
XXXS	XXS	XS	S	M	L	XS
15	15	15	15	17	19	
Ø (mm)	1,15	1,30	1,45	1,60	1,75	

Chart 2

Pathfinder

Fig. 1

Dentatus Conical Reamer

Fig. 2



Fig. 3


Diamond

Alcohol

Dry

Fig. 4


Etch

H₂O


Dry

Fig. 5


Bonding Agent

Fig. 6


Flowable Resin

Fig. 7

Post Cementation

Fig. 8

Final Post & Core Restoration

Fig. 9


Fig. 10

	CE marking		Manufacturer		Consult instructions for use		Date of manufacture
	Ultrasonic bath		Lot number		Washer-disinfector for thermal disinfection		Catalog number
	Unique device identification		Do not re-use		Medical device		Medical prescription only
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified		Caution				

SSCP is available in the European database on medical devices (EUDAMED), website: <https://ec.europa.eu/tools/eudamed> by the Basic UDI-DI for **Luscent Starter Kit**: 735008149LUC-K395 and the Basic UDI-DI for **Luscent Refills**: 735008149LUCJL.

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