

DENTATUS POST SYSTEMS - OVERVIEW

Dentatus Classic Surtex® Posts

The Dentatus Classic Surtex post has a micro-texturized metal finish, creating two-fold increase of retention compared to machined surfaces. The treated metal surface increases retention to acrylic resins and significantly reduces micro leakage. Their matte finish reduces glare for better aesthetics.

Corresponding Reamers

Helix Classic: RSA #1-6, RSB #1-6, RSC #1-5
Standard Classic: RUA #1-6, RUB #1-6

TITANIUM, SIZE #1-6

GOLD PLATED, SIZE #1-6

STAINLESS STEEL, SIZE #1-6



Dentatus Luscent Anchors®

Luscent Anchors non-metallic posts transmit light for curing composites within the deep confines of canals providing anchorage and flawless aesthetics. Resin reinforced fiberglass Luscent Anchors fulfill dentistry's preference for very strong, non-metallic posts that are consistent with conventional adhesive techniques and materials to support core foundations and aesthetic restorations. The radioluscent anchors transmit tooth color without being visible under the most adverse light conditions. They also provide a sound core-crown gingival interface radiating natural color through soft and hard tissues surrounding teeth.

Corresponding Reamers

Helix Conical: RSG #3-6, RSH #3-6, RSI #3-5
Standard Conical: RUG #1-6, RUH #4-6

XXXS: LENGTH 15 MM & Ø 1,00 MM


XXS: LENGTH 15 MM & Ø 1,15 MM

XS: LENGTH 15 MM & Ø 1,30 MM

S: LENGTH 15 MM & Ø 1,45 MM

M: LENGTH 17 MM & Ø 1,60 MM

L: LENGTH 19 MM & Ø 1,75 MM



Dentatus Twin Luscent Anchors®

Twin Luscent Anchors offer additional advantages to those of Luscent Anchors – improved mechanical retention, both against pulling and rotational forces. It gives a possibility to choose between tapered or parallel post shape, depending on which end is inserted in apical direction.

Corresponding Reamers

Helix Classic: RSA #3-6, RSB #3-6, RSC #3-5
Standard Classic: RUA #3-6, RUB #3-6

XS: LENGTH 16,7 MM & Ø 1,26 MM

S: LENGTH 16,7 MM & Ø 1,40 MM

M: LENGTH 16,7 MM & Ø 1,54 MM

L: LENGTH 16,7 MM & Ø 1,68 MM



Dentatus Luminex® Posts

Luminex light-transmitting plastic post polymerizes composites inside weak, debilitated root structures leaving an ideal space for subsequent placement of prefabricated or custom cast posts and cores. Either one of the Luminex Smooth and Grooved Plastic Posts can be used as burnout posts for custom cast posts and cores.

Corresponding Reamers

Helix Classic: RSA #1-6, RSB #1-6, RSC #1-5
Standard Classic: RUA #1-6, RUB #1-6





SMOOTH, SIZE #1-6

GROOVED, SIZE #1-6



CORRESPONDING REAMERS

READ MORE ON PAGE 20

Post system & Post Sizes		Corresponding Reamers								
Surtex Luminex	Twin Luscent		Helix Classic Reamers				Standard Classic Reamers			
			Short 28 mm	Long 33 mm	X Long 37 mm		Short 28 mm	Long 33 mm		
			#1	-	RSA-1/3		RSB-1/3	RSC-1/2	RUA-1/3	RUB-1/3
			#2	-	RSA-2/3		RSB-2/3	RSC-2/2	RUA-2/3	RUB-2/3
			#3	White (XS)	RSA-3/3		RSB-3/3	RSC-3/2	RUA-3/3	RUB-3/3
			#4	Yellow (S)	RSA-4/3		RSB-4/3	RSC-4/2	RUA-4/3	RUB-4/3
			#5	Red (M)	RSA-5/3		RSB-5/3	RSC-5/2	RUA-5/3	RUB-5/3
#6	Blue (L)	RSA-6/3	RSB-6/3	-	RUA-6/3	RUB-6/3				
	Luscent		Helix Conical Reamers				Standard Conical Reamers			
	Purple (XXXS)		-	-	-		RUG-1/2			
	Green (XXS)		-	-	-		RUG-2/2			
	White (XS)		RSG-3/3	RSH-3/3	RSI-3/2		RUG-3/2			
	Yellow (S)		RSG-4/3	RSH-4/3	RSI-4/2		RUG-4/2	RUH-4/2		
	Red (M)		RSG-5/3	RSH-5/3	RSI-5/2		RUG-5/2	RUH-5/2		
	Blue (L)		RSG-6/3	RSH-6/3	-		RUG-6/2	RUH-6/2		

DENTATUS CLASSIC SURTEX® POSTS

PREMIUM DENTAL POSTS WITH PROVEN RELIABILITY

Superior Retention Through Design

The Advantages

- Surtex® surface treatment
- Superior retention
- Assortment box for ease-of-use
- Manufactured in Sweden
- 100% quality control


The original Dentatus Classic Surtex® Posts have become a global standard in dental anchors, and offer dentists a reliable and simple solution for post and core build-up. Our high-quality, pre-fabricated posts are dependable and economic for most endodontic and prosthodontic treatment needs, available in a range of materials including pure titanium, stainless steel and gold-plated.

Whatever your crown retention requirements, Dentatus Classic Surtex® Posts has the solution. Available in titanium, stainless steel and goldplated, together with our dedicated reamers they form a state-of-the-art system for post retention in endodontically treated root canals. Dentatus Classic Surtex Posts – a continuously developed global standard since 1932.

Square crosscut head for core material adaption

Parallell sides and surface treated for extra retention

Anatomical geometry for tapered roots



DENTATUS CLASSIC SURTEX® POST SIZES

	SHORT (7,8 MM)	MEDIUM (9,3 MM)	LONG (11,8 MM)	EXTRA LONG (14,2 MM)	SUPER LONG (17 MM)
1	Ø 1,05 mm	Ø 1,05 mm	Ø 1,05 mm		
2	Ø 1,20 mm	Ø 1,20 mm	Ø 1,20 mm		
3	Ø 1,35 mm	Ø 1,35 mm	Ø 1,35 mm	Ø 1,35 mm	Ø 1,35 mm
4	Ø 1,50 mm	Ø 1,50 mm	Ø 1,50 mm	Ø 1,50 mm	Ø 1,50 mm
5	Ø 1,65 mm	Ø 1,65 mm	Ø 1,65 mm	Ø 1,65 mm	Ø 1,65 mm
6	Ø 1,80 mm	Ø 1,80 mm	Ø 1,80 mm	Ø 1,80 mm	Ø 1,80 mm



Available in Titanium, Gold Plated and Stainless Steel

STARTER KITS

CLASSIC SURTEX TITANIUM

Assortment Kit

- TST-60** 60 Assorted Posts (20 sizes, 3 posts of each)
 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
 1 Hollow Key
 1 Cross Key
 20 Stoppers



Complete Starter Kit

- TST-B** 120 Assorted Posts (20 sizes, 6 posts of each)
 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
 1 Hollow Key
 1 Cross Key
 20 Stoppers



Starter Refill Kit - No Reamers

- TST** 120 Assorted Posts (20 sizes, 6 posts of each)
 1 Hollow Key
 1 Cross Key
 20 Stoppers



STARTER KITS

CLASSIC SURTEX GOLD PLATED

Assortment Kit

- RST-60** 60 Assorted Posts (20 sizes, 3 posts of each)
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers



Complete Starter Kit

- RST-B** 120 Assorted Posts (20 sizes, 6 posts of each)
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers



Starter Refill Kit - No Reamers

- RST** 120 Assorted Posts (20 sizes, 6 posts of each)
1 Hollow Key
1 Cross Key
20 Stoppers



CLASSIC SURTEX STAINLESS STEEL

Assortment Kit

- SST-60** 60 Assorted Posts (20 sizes, 3 posts of each)
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers



Complete Starter Kit

- SST-B** 120 Assorted Posts (20 sizes, 6 posts of each)
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers



Starter Refill Kit - No Reamers

- SST** 120 Assorted Posts (20 sizes, 6 posts of each)
1 Hollow Key
1 Cross Key
20 Stoppers



REFILLS

CLASSIC SURTEX TITANIUM

Short, 7.8 mm, 15 pcs/pkg

TST-S1	Diameter 1.05 mm
TST-S2	Diameter 1.20 mm
TST-S3	Diameter 1.35 mm
TST-S4	Diameter 1.50 mm
TST-S5	Diameter 1.65 mm
TST-S6	Diameter 1.80 mm



Medium, 9.3 mm, 15 pcs/pkg

TST-M1	Diameter 1.05 mm
TST-M2	Diameter 1.20 mm
TST-M3	Diameter 1.35 mm
TST-M4	Diameter 1.50 mm
TST-M5	Diameter 1.65 mm
TST-M6	Diameter 1.80 mm



Long, 11.8 mm, 15 pcs/pkg

TST-L1	Diameter 1.05 mm
TST-L2	Diameter 1.20 mm
TST-L3	Diameter 1.35 mm
TST-L4	Diameter 1.50 mm
TST-L5	Diameter 1.65 mm
TST-L6	Diameter 1.80 mm



Extra Long, 14.2 mm, 15 pcs/pkg

TST-EL3	Diameter 1.35 mm
TST-EL4	Diameter 1.50 mm
TST-EL5	Diameter 1.65 mm
TST-EL6	Diameter 1.80 mm



Super Long, 17 mm, 15 pcs/pkg

TST-SL3	Diameter 1.35 mm
TST-SL4	Diameter 1.50 mm
TST-SL5	Diameter 1.65 mm
TST-SL6	Diameter 1.80 mm



REFILLS

CLASSIC SURTEX GOLD PLATED

Short, 7.8 mm, 15 pcs/pkg

RST-S1	Diameter 1.05 mm
RST-S2	Diameter 1.20 mm
RST-S3	Diameter 1.35 mm
RST-S4	Diameter 1.50 mm
RST-S5	Diameter 1.65 mm
RST-S6	Diameter 1.80 mm



Medium, 9.3 mm, 15 pcs/pkg

RST-M1	Diameter 1.05 mm
RST-M2	Diameter 1.20 mm
RST-M3	Diameter 1.35 mm
RST-M4	Diameter 1.50 mm
RST-M5	Diameter 1.65 mm
RST-M6	Diameter 1.80 mm



Long, 11.8 mm, 15 pcs/pkg

RST-L1	Diameter 1.05 mm
RST-L2	Diameter 1.20 mm
RST-L3	Diameter 1.35 mm
RST-L4	Diameter 1.50 mm
RST-L5	Diameter 1.65 mm
RST-L6	Diameter 1.80 mm



Extra Long, 14.2 mm, 15 pcs/pkg

RST-EL3	Diameter 1.35 mm
RST-EL4	Diameter 1.50 mm
RST-EL5	Diameter 1.65 mm
RST-EL6	Diameter 1.80 mm



Super Long, 17 mm, 15 pcs/pkg

RST-SL3	Diameter 1.35 mm
RST-SL4	Diameter 1.50 mm
RST-SL5	Diameter 1.65 mm
RST-SL6	Diameter 1.80 mm



REFILLS

CLASSIC SURTEX STAINLESS STEEL

Short, 7.8 mm, 12 pcs/pkg

SST-S1	Diameter 1.05 mm
SST-S2	Diameter 1.20 mm
SST-S3	Diameter 1.35 mm
SST-S4	Diameter 1.50 mm
SST-S5	Diameter 1.65 mm
SST-S6	Diameter 1.80 mm



Medium, 9.3 mm, 12 pcs/pkg

SST-M1	Diameter 1.05 mm
SST-M2	Diameter 1.20 mm
SST-M3	Diameter 1.35 mm
SST-M4	Diameter 1.50 mm
SST-M5	Diameter 1.65 mm
SST-M6	Diameter 1.80 mm



Long, 11.8 mm, 12 pcs/pkg

SST-L1	Diameter 1.05 mm
SST-L2	Diameter 1.20 mm
SST-L3	Diameter 1.35 mm
SST-L4	Diameter 1.50 mm
SST-L5	Diameter 1.65 mm
SST-L6	Diameter 1.80 mm



Extra Long, 14.2 mm, 12 pcs/pkg

SST-EL3	Diameter 1.35 mm
SST-EL4	Diameter 1.50 mm
SST-EL5	Diameter 1.65 mm
SST-EL6	Diameter 1.80 mm



Super Long, 17 mm, 12 pcs/pkg

SST-SL3	Diameter 1.35 mm
SST-SL4	Diameter 1.50 mm
SST-SL5	Diameter 1.65 mm
SST-SL6	Diameter 1.80 mm



ACCESSORIES

FOR DENTATUS METAL POSTS

NLH 1 Hollow Key

NLK 1 Cross Key

