

# Dentatus Reamers

Dentatus

CE  
0413

Dentatus AB is certified to ISO 13485

 **Dentatus AB**  
Finspångsgatan 42  
SE-163 53 Spånga, Sweden



**SE** +46 8 546 509 00  
info@dentatus.se

**US** + 1 212 481 1010  
dentatus@dentatus.com

[www.dentatus.com](http://www.dentatus.com)

MADE IN SWEDEN

BL-6611-03\_2019-11-22

## Dentatus Reamers



**Dentatus Helix Reamers**  
With sharp cutting edges, they have high cutting capacity, smooth cutting action and improved transport of cutting debris. Available in Classic and Conical shapes.

**Dentatus Standard Reamers**  
The Original Dentatus reamer. Available in Classic and Conical shapes.

**Dentatus Special Reamers**  
The Probos I Pathfinder is ideal for initial and safe removal of root filling material. The tip design minimizes the risk for perforation of the dental walls.

The Probos II Reamer/Router increases the surface by indenting the dental wall to improve retention, when used laterally in the canal.

Read more on [www.dentatus.com](http://www.dentatus.com)

# POSTS & CORRESPONDING REAMERS

Post system (post sizes)

Corresponding Reamers

For assortments - contact your local dealer

Surtex  
Luminex

Twin Luscent

**Helix**



Classic Reamers	
Short 28 mm	Long 33 mm
X Long 37 mm	
RSA-1/3	RSB-1/3
RSA-2/3	RSB-2/3
RSA-3/3	RSB-3/3
RSA-4/3	RSB-4/3
RSA-5/3	RSB-5/3
RSA-6/3	RSB-6/3
RSC-1/2	RSC-2/2
RSC-3/2	RSC-4/2
RSC-5/2	



**Standard**

Classic Reamers

Classic Reamers	
Short 28 mm	Long 33 mm
RUA-1/3	RUB-1/3
RUA-2/3	RUB-2/3
RUA-3/3	RUB-3/3
RUA-4/3	RUB-4/3
RUA-5/3	RUB-5/3
RUA-6/3	RUB-6/3

- #1
- #2
- #3
- #4
- #5
- #6

- 
- 
- White (XS)
- Yellow (S)
- Red (M)
- Blue (L)

Luscent

Conical Reamers



- Purple (XXXS)
- Green (XXS)
- White (XS)
- Yellow (S)
- Red (M)
- Blue (L)

-	-	-
-	-	-
RSG-3/3	RSH-3/3	RSI-3/2
RSG-4/3	RSH-4/3	RSI-4/2
RSG-5/3	RSH-5/3	RSI-5/2
RSG-6/3	RSH-6/3	-

Conical Reamers

- RUG-1/2
- RUG-2/2
- RUG-3/2
- RUG-4/2
- RUG-5/2
- RUG-6/2
- RUH-4/2
- RUH-5/2
- RUH-6/2