

DENTATUS AB

A world leader in the design & manufacture of dental products



Dentatus Profin[®] Directional System

The ultimate instrumentation for finishing and polishing for professionals

- Contra-angle Profin PDX
- Hand instrument PDH
- Lamineer Tips
- Per-lo-Tor Instruments
- Made in Sweden



Dentatus Profin[®]

Your Companion to Dental Perfection

Perfect for:

- Hygienists
- Orthodontists
- Periodontists
- Prosthodontists
- Implantologists



Dentatus Profin[®] PDX (1/2)

- 1.2 mm reciprocating tip movement
- Fixed or rotating mode
- For hard-to-reach areas
- No heat build-up
- Can be used close to the gingiva




Dentatus Profin[®] PDX (2/2)


- Tactile ergonomics
- Smooth glare free surface
- Gentle cleaning of root surfaces (with Per-lo-Tor tips)
- Made in Sweden




The unique reciprocating movement by Profin[®] PDX



Rotating instrumentation



Reciprocating instrumentation



The Profin movement


Profin[®] for Hygienists & Periodontists

- Gentle cleaning of root surfaces with minimal or no damage to the root cementum (with Per-lo-Tor tips)
- Efficient removal of restorational overhangs
- Professional cleaning and polishing of all dental surfaces



Profin[®] for Orthodontists

- Interproximal reduction during orthodontic treatment with fixed appliances
- Clean-up of tooth surfaces after de-bonding
- Stripping with IPR tips



Profin[®] for Prosthodontics





- Fine-tuning of all kinds of preparations
- Individually adapted shaping and contouring for enhanced aesthetics
- Finishing and polishing of all restorations, also in hard-to-reach areas
- Wide area of LTA Lamineer[®] & EVA tips



Profin[®] for Implantologists

- Single use plastic EVA tips
- Use with dedicated cleaning past
- Gentle to prosthetics and gingiva
- No heat build-up
- Gentle to implants








Profin[®] PDH

Achieve the Highest Levels of Dental Finishing

- Based on the established Profin/EVA system
- The same wide range of interchangeable diamond Lamineer tips
- For safe and precise finishing work
- For posterior and anterior usage



Dentatus Profin[®] PDH




Dentatus LTA Lamineer[®] Tips

Reciprocating tools for finishing and polishing of all tooth surfaces



LTA



IPR



Dentatus LTA Lamineer[®] Tips

- Provides excellent and well defined surfaces, also in hard-to-reach areas
- Color coded abrasive diamond tips in different shapes, grits and sizes for specific all kinds of treatments
- Safe-sided to prevent damage to adjacent tissue or enamel
- Sterilizable to provide long term economical use

Dentatus LTA Lamineer[®] Tips

Additional Shapes and Functions:

- Shorter tips for CAD/CAM restorations
- EVA plastic tips for Implant Cleaning
- IPR tips for Interproximal reduction
- Longer, more flexible tips for contours
- Profin/EVA tips can also be used with W&H and NSK contra-angles for Profin
- ...and more

Dentatus Per-lo-Tor™ Instruments

The ultimate tools for safe and mechanical root-planing with Profin PDX



Per-lo-Tor






Prevents Further Removal of Root Cementum Once Surface is Clean


- For mechanical tooth cleaning
- Safe root-planing with Profin
- Developed by Prof. Per Axelsson
- Made of high quality stainless steel
- Made for long lasting performance



Why use Per-lo-Tor



The curette



Ultrasonic/
High Speed Scaler



Consequence of recurrence of
plaque on partially exposed root
dentin.

DENTATUS AB

Thank you for your attention.