

DENTATUS PROFIN IPR TIPS FOR PROFIN/EVA

FOR SAFE INTERPROXIMAL REDUCTION

Increase the Durability of Orthodontic Treatment Results

The Advantages

- Safety is top concern for the patients
- Control is in your hand
- Performance made easier
- Esthetics is all natural

The Profin IPR™ tips for interproximal reduction are developed to use together with the highly acclaimed reciprocating Profin Directional System.


Profin is an essential addition to any discipline in dentistry and, now more than ever, in orthodontics. For a dentist, patient safety is top priority; that is why we developed the Profin IPR. With Profin IPR, dentists can correct disproportionate tooth sizes, eliminate crowding, protect enamel and increase the durability of orthodontic treatment results. Profin IPR diminishes the risk of excessive stripping thereby ensuring more conservative removal of enamel in prescribed amounts.

Profin IPR also eliminates the risk of damage to gingival tissue.

Safe sided tips

Double sided tips

Color coded tips



Safety is Top Concern for Our Patients

Profin IPR diminishes the risk of excessive stripping thereby ensuring more conservative removal of enamel in prescribed amounts. Profin IPR also eliminates the risk of damage to gingival tissue.

Control is in Your Hands

Used with the reciprocating Profin Directional System, Profin IPR's movement is fluid as the 1.2 mm axial stroke allows for concise corrections and safe, quick and efficient modification


Performance Made Easier

Profin IPR's controlled reduction avoids rounding teeth. It allows you to reach all the right places for optimal results, shaping, finishing and preserving the natural anatomy of the tooth.



Esthetics Have Never Been More Natural

With Profin IPR practitioners can easily recreate teeth's natural contours with a variety of Lamineer Tips in numerous shapes, sizes and grits.



Product discription:

LTA-P25R: Violet shank, safe-sided – abrasive on right side when shank points down and back of tip towards user. Blade thickness 0.25mm.

LTA-P25L: Violet shank, safe-sided – abrasive on left side when shank points down and back of tip towards user. Blade thickness 0.25mm.

LTA-P3: Yellow shank. Double sided. Blade thickness 0.3mm.

LTA-P4: Green shank. Double sided. Blade thickness 0.4mm.

LTA-P5: Red shank. Double sided. Blade thickness 0.5mm.

All IPR tips have the same grit – equivalent to 50 μ . Color coding is relating to blade thickness only.

Starter kit and refills

You can buy the Profin IPR kit with 5 tips and the stand or you can buy the tips separately in packs of 2 pcs.



Dentatus AB
www.dentatus.com
SE: +46 8 546 509 00
info@dentatus.se
US: +1 212 481 1010
dentatus@dentatus.com

