

DENTATUS EVA TIPS

FOR IMPLANT CLEANING

Efficient and Gentle Cleaning

The Advantages


- Single use plastic tips
- Use with dedicated cleaning paste
- No heat build-up
- Gentle to implants
- Gentle to prosthetics and gingiva

By using the plastic EVA tips for cleaning implants and implant retained prosthetics, an efficient and gentle cleaning can be provided. The soft plastic tips do not damage the implant surface and are gentle both against the implant, the prosthesis and against the patient's gums. The gentle reciprocating movement of 1.2 mm allows for easy and efficient cleaning without heat development. It makes it easier for the dentist or dental hygienist to access hard to reach areas and provide a more comfortable experience for the patient.


The plastic tips are disposable and used together with appropriate polishing paste.



Profin PDX hand piece




EVA 5000 plastic tip



Profin PDH hand instrument

The green tip: EVA 2000 has a sharp pointed tip for cleaning of all approximal areas.

The blue tip: EVA 5000, thin spatula with round tip.



Dentatus AB
www.dentatus.com
SE: +46 8 546 509 00
info@dentatus.se
US: +1 212 481 1010
dentatus@dentatus.com



BL6712EN-03_2018-10