

DENTATUS LUMINEX® POST SYSTEM

FOR RESTORING THIN-WALLED ROOTS

The Visible Cure for Saving Compromised Roots

The Advantages

- Light guides for building up roots
- To strengthen debilitated roots
- Create a “standard” root canal
- Easy single office visit procedure

The Luminex® post allows the polymerizing light to pass down through light cured composites inside weak, debilitated root structures, in order to strengthen the root and to leave an ideal space for subsequent placement of prefabricated or custom cast posts and cores.

Markings for different sizes

Grooved or smooth surface

Light transmitting



REFILLS

LUMINEX SMOOTH POSTS

Long, 22 mm, 20 pcs/pkg

DIP-1	Diameter 1.05 mm
DIP-2	Diameter 1.20 mm
DIP-3	Diameter 1.35 mm
DIP-4	Diameter 1.50 mm
DIP-5	Diameter 1.65 mm
DIP-6	Diameter 1.80 mm



LUMINEX GROOVED BURNOUT POSTS

Long, 22 mm, 20 pcs/pkg

DCP-1	Diameter 1.05 mm
DCP-2	Diameter 1.20 mm
DCP-3	Diameter 1.35 mm
DCP-4	Diameter 1.50 mm
DCP-5	Diameter 1.65 mm
DCP-6	Diameter 1.80 mm

