

DENTATUS PER-IO-TOR™ INSTRUMENTS

THE ULTIMATE TOOLS FOR SAFE MECHANICAL
ROOT-PLANING

Prevents Further Removal of Root Cementum Once Surface is Clean


The Advantages

- For mechanical tooth cleaning
- Safe root-planing with Profin
- Developed by Prof. Per Axelsson
- Made of high quality stainless steel
- Made for long lasting performance

Used in the reciprocating Profin handpieces, Per-Io-Tor tips gently smoothen the root surface without over-instrumentation.

For Profin reciprocating hand piece

All instruments are made of high quality stainless steel



FACTS PER-IO-TOR FOR ROOT PLANING:

- All Per-lo-Tor instruments are made of high quality stainless steel and may be cleaned in ultrasonic baths, in thermal washers / disinfectors and in autoclaves at 134° C.
- Per-lo-Tor 1-4 may be used in all Profin and EVA handpieces.
- Per-lo-Tor 5-6 should be used in a Profin handpiece, since these instruments are used in a locked position.
- Degree of removal: approximately 0,005 mm at 2-3 minutes use until the root cementum is smooth.
- All Per-lo-Tor instruments should be used at approximately 10.000 rpm.
- The Per-lo-Tor instruments have been developed by Prof. Per Axelsson, Sweden, and are carefully documented.

STARTER KIT

PER-IO-TOR, RECIPROCATING

Instruments and Stand

TOR-K Tip stand with 6 Per-lo-Tor instruments.
One of each of TOR #1 to #6.



REFILLS

PER-IO-TOR, RECIPROCATING

Instruments, 1pc/pkg

TOR-1	#1, PER-IO-TOR Tip, Round with grooves, small	
TOR-2	#2, PER-IO-TOR Tip, Round with grooves, medium	
TOR-3	#3, PER-IO-TOR Tip, Spatula with eyelet	
TOR-4	#4, PER-IO-TOR Tip, Spatula with ridges	
TOR-5	#5, PER-IO-TOR Tip, Concave	
TOR-6	#6, PER-IO-TOR Tip, Convex	