

DENTATUS POST SYSTEMS - OVERVIEW

Dentatus Classic Surtex® Posts

The Dentatus Classic Surtex post has a micro-texturized metal finish, creating two-fold increase of retention compared to machined surfaces. The treated metal surface increases retention to acrylic resins and significantly reduces micro leakage. Their matte finish reduces glare for better aesthetics.

Corresponding Reamers

Helix Classic: RSA #1-6, RSB #1-6, RSC #1-5

Standard Classic: RUA #1-6, RUB #1-6

TITANIUM, SIZE #1-6

GOLD PLATED, SIZE #1-6

STAINLESS STEEL, SIZE #1-6



Dentatus Luscent Anchors®

Luscent Anchors non-metallic posts transmit light for curing composites within the deep confines of canals providing anchorage and flawless aesthetics. Resin reinforced fiberglass Luscent Anchors fulfill dentistry's preference for very strong, non-metallic posts that are consistent with conventional adhesive techniques and materials to support core foundations and aesthetic restorations. The radioluscent anchors transmit tooth color without being visible under the most adverse light conditions. They also provide a sound core-crown gingival interface radiating natural color through soft and hard tissues surrounding teeth.

Corresponding Reamers

Helix Conical: RSG #3-6, RSH #3-6, RSI #3-5

Standard Conical: RUG #1-6, RUH #4-6

XXS: LENGTH 15 MM & Ø 1,00 MM

XXS: LENGTH 15 MM & Ø 1,15 MM

XS: LENGTH 15 MM & Ø 1,30 MM

S: LENGTH 15 MM & Ø 1,45 MM

M: LENGTH 17 MM & Ø 1,60 MM

L: LENGTH 19 MM & Ø 1,75 MM



Dentatus Twin Luscent Anchors®

Twin Luscent Anchors offer additional advantages to those of Luscent Anchors – improved mechanical retention, both against pulling and rotational forces. It gives a possibility to choose between tapered or parallel post shape, depending on which end is inserted in apical direction.

Corresponding Reamers

Helix Classic: RSA #3-6, RSB #3-6, RSC #3-5


Standard Classic: RUA #3-6, RUB #3-6

XS: LENGTH 16,7 MM & Ø 1,26 MM

S: LENGTH 16,7 MM & Ø 1,40 MM

M: LENGTH 16,7 MM & Ø 1,54 MM

L: LENGTH 16,7 MM & Ø 1,68 MM



Dentatus Luminex® Posts

Luminex light-transmitting plastic post polymerizes composites inside weak, debilitated root structures leaving an ideal space for subsequent placement of prefabricated or custom cast posts and cores. Either one of the Luminex Smooth and Grooved Plastic Posts can be used as burnout posts for custom cast posts and cores.

Corresponding Reamers

Helix Classic: RSA #1-6, RSB #1-6, RSC #1-5

Standard Classic: RUA #1-6, RUB #1-6





SMOOTH, SIZE #1-6

GROVED, SIZE #1-6



CORRESPONDING REAMERS

READ MORE ON PAGE 20

Post system & Post Sizes		Corresponding Reamers									
Surtex Luminex	Twin Luscent		Helix Classic Reamers				Standard Classic Reamers				
			Short 28 mm	Long 33 mm	X Long 37 mm		Short 28 mm	Long 33 mm			
	#1		-	RSA-1/3	RSB-1/3		RSC-1/2	RUA-1/3	RUB-1/3		
	#2		-	RSA-2/3	RSB-2/3		RSC-2/2	RUA-2/3	RUB-2/3		
	#3		White (XS)	RSA-3/3	RSB-3/3		RSC-3/2	RUA-3/3	RUB-3/3		
	#4		Yellow (S)	RSA-4/3	RSB-4/3		RSC-4/2	RUA-4/3	RUB-4/3		
	#5		Red (M)	RSA-5/3	RSB-5/3		RSC-5/2	RUA-5/3	RUB-5/3		
	#6		Blue (L)	RSA-6/3	RSB-6/3		-	RUA-6/3	RUB-6/3		
			Luscent		Helix Conical Reamers				Standard Conical Reamers		
	Purple (XXXS)	-	-		-	RUG-1/2					
	Green (XXS)	-	-		-	RUG-2/2					
	White (XS)	RSG-3/3	RSH-3/3		RSI-3/2	RUG-3/2					
	Yellow (S)	RSG-4/3	RSH-4/3		RSI-4/2	RUG-4/2	RUH-4/2				
	Red (M)	RSG-5/3	RSH-5/3		RSI-5/2	RUG-5/2	RUH-5/2				
	Blue (L)	RSG-6/3	RSH-6/3	-	RUG-6/2	RUH-6/2					

DENTATUS REAMERS

CREATING THE PERFECT POSTULATE

The Perfect Match to Dentatus Post Systems

The Advantages


- For a perfect match between the root and the canal posts
- High quality stainless steel
- Wide choice of sizes

Dentatus reamers are designed to create a perfect match between dental root canals and the shapes and sizes of the corresponding posts. Made in Sweden from surgical high quality steel for high quality performance.

Reliable and long lasting performance

Sizes are indicated by numbers of grooves on the shank

Sharp cutting edges to create perfect seat for Dentatus Posts



ASSORTED PACKAGES

HELIX CLASSIC REAMERS

Assorted Packages, 3 pcs/pkg

RSA-123	Helix Classic Reamers, Short (1 of each size 1, 2 ,3)
RSA-345	Helix Classic Reamers, Short (1 of each size 3, 4 ,5)
RSA-456	Helix Classic Reamers, Short (1 of each size 4, 5, 6)
RSB-123	Helix Classic Reamers, Long (1 of each size 1, 2 ,3)
RSB-345	Helix Classic Reamers, Long (1 of each size 3, 4 ,5)
RSB-456	Helix Classic Reamers, Long (1 of each size 4, 5, 6)
RSC-123	Helix Classic Reamers, Extra Long (1 of each size 1, 2 ,3)
RSC-345	Helix Classic Reamers, Extra Long (1 of each size 3, 4 ,5)



REFILLS

HELIX CLASSIC REAMERS

Short, 28 mm, 3 pcs/pkg

RSA-1/3	(#1) XXX-Small	Diameter 1.05 mm
RSA-2/3	(#2) XX-Small	Diameter 1.20 mm
RSA-3/3	(#3) X-Small	Diameter 1.35 mm
RSA-4/3	(#4) Small	Diameter 1.50 mm
RSA-5/3	(#5) Medium	Diameter 1.65 mm
RSA-6/3	(#6) Large	Diameter 1.80 mm



Long, 33 mm, 3 pcs/pkg

RSB-1/3	(#1) XXX-Small	Diameter 1.05 mm
RSB-2/3	(#2) XX-Small	Diameter 1.20 mm
RSB-3/3	(#3) X-Small	Diameter 1.35 mm
RSB-4/3	(#4) Small	Diameter 1.50 mm
RSB-5/3	(#5) Medium	Diameter 1.65 mm
RSB-6/3	(#6) Large	Diameter 1.80 mm



Extra Long, 37 mm, 2 pcs/pkg

RSC-1/2	(#1) XXX-Small	Diameter 1.05 mm
RSC-2/2	(#2) XX-Small	Diameter 1.20 mm
RSC-3/2	(#3) X-Small	Diameter 1.35 mm
RSC-4/2	(#4) Small	Diameter 1.50 mm
RSC-5/2	(#5) Medium	Diameter 1.65 mm



ASSORTED PACKAGES

HELIX CONICAL REAMERS

Assorted Packages, 3 pcs/pkg

RSG-345	Helix Conical Reamers, Short (1 of each size 3, 4 ,5)
RSG-456	Helix Conical Reamers, Short (1 of each size 4, 5, 6)
RSH-345	Helix Conical Reamers, Long (1 of each size 3, 4 ,5)
RSH-456	Helix Conical Reamers, Long (1 of each size 4, 5, 6)
RSI-345	Helix Conical Reamers, Extra Long (1 of each size 3, 4 ,5)



REFILLS

HELIX CONICAL REAMERS

Short, 28 mm, 3 pcs/pkg

RSG-3/3	(#3) X-Small	Diameter 1.35 mm
RSG-4/3	(#4) Small	Diameter 1.50 mm
RSG-5/3	(#5) Medium	Diameter 1.65 mm
RSG-6/3	(#6) Large	Diameter 1.80 mm



Long, 33 mm, 3 pcs/pkg

RSH-3/3	(#3) X-Small	Diameter 1.35 mm
RSH-4/3	(#4) Small	Diameter 1.50 mm
RSH-5/3	(#5) Medium	Diameter 1.65 mm
RSH-6/3	(#6) Large	Diameter 1.80 mm



Extra Long, 37 mm, 2 pcs/pkg

RSI-3/2	(#3) X-Small	Diameter 1.35 mm
RSI-4/2	(#4) Small	Diameter 1.50 mm
RSI-5/2	(#5) Medium	Diameter 1.65 mm



ASSORTED PACKAGES

STANDARD CLASSIC REAMERS

Assorted Packages, 6 pcs/pkg

RUA-1/6	Standard Classic Reamers, Short (1 of each size)
RUB-1/6	Standard Classic Reamers, Long (1 of each size)



REFILLS

STANDARD CLASSIC REAMERS

Short, 28 mm, 3 pcs/pkg

RUA-1/3FP	(#1) XXX-Small	Diameter 1.05 mm
RUA-2/3FP	(#2) XX-Small	Diameter 1.20 mm
RUA-3/3FP	(#3) X-Small	Diameter 1.35 mm
RUA-4/3FP	(#4) Small	Diameter 1.50 mm
RUA-5/3FP	(#5) Medium	Diameter 1.65 mm
RUA-6/3FP	(#6) Large	Diameter 1.80 mm



Long, 33 mm, 3 pcs/pkg

RUB-1/3FP	(#1) XXX-Small	Diameter 1.05 mm
RUB-2/3FP	(#2) XX-Small	Diameter 1.20 mm
RUB-3/3FP	(#3) X-Small	Diameter 1.35 mm
RUB-4/3FP	(#4) Small	Diameter 1.50 mm
RUB-5/3FP	(#5) Medium	Diameter 1.65 mm
RUB-6/3FP	(#6) Large	Diameter 1.80 mm



STANDARD CONICAL REAMERS

Short, 28 mm, 2 pcs/pkg

RUG-1/2	(#1) XXX-Small	Diameter 1.05 mm
RUG-2/2	(#2) XX-Small	Diameter 1.20 mm
RUG-3/2	(#3) X-Small	Diameter 1.35 mm
RUG-4/2	(#4) Small	Diameter 1.50 mm
RUG-5/2	(#5) Medium	Diameter 1.65 mm
RUG-6/2	(#6) Large	Diameter 1.80 mm



Long, 33 mm, 2 pcs/pkg

RUH-4/2	(#4) Small	Diameter 1.50 mm
RUH-5/2	(#5) Medium	Diameter 1.65 mm
RUH-6/2	(#6) Large	Diameter 1.80 mm



ASSORTED PACKAGES

SPECIAL REAMERS

Assorted Packages

RUF-1/6	6 Probos® I, Pathfinder Assorted (1 of each size)
RUD	4 Probos® II, Reamer/Router Assorted (1 of each size)

REFILLS

PROBOS I, PATHFINDER

Long, 33 mm, 3 pcs/pkg

RUF-1/3	#1	Diameter 0.60 mm
RUF-2/3	#2	Diameter 0.75 mm
RUF-3/3	#3	Diameter 0.90 mm
RUF-4/3	#4	Diameter 1.05 mm
RUF-5/3	#5	Diameter 1.20 mm
RUF-6/3	#6	Diameter 1.35 mm



PROBOS II, REAMER/ROUTER

Short, 28 mm, 3 pcs/pkg

RUD-3/3	#3	Diameter 1.35 mm
RUD-4/3	#4	Diameter 1.50 mm
RUD-5/3	#5	Diameter 1.65 mm
RUD-6/3	#6	Diameter 1.80 mm



MANDRELS FOR CUTTING AND POLISHING DISCS

Well balanced high quality mandrels for straight (HP) and angled hand pieces (RA) have been manufactured by Dentatus for many years. Used together with cutting and polishing discs both at the surgery and at the dental lab.



MANDRELS

Diameter 2,35mm, 12pcs/pkg

MLR-RA	Dentatus Mandrels RA, Length 24 mm Diameter 2,35mm
MLR-HP	Dentatus Mandrels, HP, Length 47 mm Diameter 2,35mm

