



In First Impressions, George Freedman, DDS, gives readers a brief summary of products that have recently been introduced to dentistry, based on his clinical experience.



ANEW Narrow Body Implants

tatus

face the following 2 problems with implant delivery: insufficient bone width and limited mesio-distal space between existing teeth. Traditional implants require a minimum of buccal and lingual support, and the tipping in of the implant around an earlier extraction site leaves less interdental space between the diameter of the implant. Dentatus' FDA-approved ANEW Narrow Body Implants System offers a fixed and removable prosthetic

options for even the most compromised cases. Designed to complement mainstream implant systems, the ANEW technique includes various implant platforms for both immediate and long-term use of screw-retained single- or multiunit prosthetics as well as for temporary restorations that offer patients a comfortable dentate state throughout their implant treatment. The Elypse Platform allows ANEW, when used with the ATLAS Denture Comfort System, to provide reliable


treatment for removable prostheses. Its retentive screw caps facilitate easier removal and reattachment of the prosthetics. ANEW's non-hygroscopic screw cap allows for implants to be monitored, altered, and adapted to a permanent restoration. It also protects grafting sites, maintains the desired soft-tissue architecture, and ensures that the final implant is located in the ideal aesthetic position. It is manufactured into a single piece of Grade-5 titanium alloy and is available in 1.8-,

2.2-, and 2.4-mm narrow-body diameters. ANEW's tapering implant improves initial stability while its square platform prevents prosthetic rotation, and its surface-treated coating promotes osseointegration, offering instant support for nonretentive failing abutments and restorations.

For more information, call us free (800) 323-3136 or visit the website dentatususa.com. ♦

FREEinfo, circle 63

When you've considered all the options, isn't it time to consider ANEW?



Dentatus
ANEW[®]
Narrow Body Implants

Dentatus ANEW implants provide gold-standard treatment options for many more patients. In areas of limited bone width, mesial-distal space or converging roots, Dentatus ANEW is often an ideal solution available in 1.8, 2.2 & 2.4mm diameters and three thread lengths. With a tapered apical end, Dentatus ANEW Implants can be placed in interdental spaces as narrow as 3.5mm without the need for bone augmentation or orthodontic interventions.

Dentatus ANEW is the only narrow-body implant with a screw-retained prosthetic system and with over 10 years of clinical research to support safe and reliable long-term use.

Expanding implantology with more treatment options.

www.dentatususa.com • 800-323-3136



©2015 Dentatus USA, Ltd. Patented and Patents Pending

KNOWN CLINICIANS TAKE THE PODIUM
Dr. Paul Petruccaro, Sept
WSP • Newport Beach
Ziv Mazor
AP • O