

Dentatus® **Matrix Bands**



Dentatus Matrix bands are manufactured from the highest quality stainless steel. The standard shape bands are available in two thicknesses: 0.050 and 0.040 mm.

The 0.040 bands are made of high quality stainless steel, and the thicker bands, 0.050, are made from stainless steel with special dead soft quality. The choice of steel quality gives the bands high strength and good adjustability. The bands provide very good prerequisites for well contoured and defined interproximal contacts.

To reduce total interproximal matrix band thickness on mesio-occlusal and disto-occlusal restorations without reducing actual band thickness, the MoDo bands are specially cut to go under the contact on the intact side of the tooth.

Dentatus quality has proven to be highly regarded by users all over the world for a long period of time.





Matrix Bands

ORIGINAL DENTATUS MATRIX BANDS

Stainless steel, Thin 0.040 mm, 12 pcs.

Stainless dead soft steel, 0.050 mm
 Package available with 12 (/DS) bands or
 100 (/100FP) bands, e.g. MDD-1/DS or MDD-1/100FP.

MDD-1RF

MDD-2RF

MDD-1

MDD-2

MDD-3RF

MDD-4RF

MDD-3







MDD-4

MDD-5RF

MDD-6RF

MDD-5

MDD-6



DENTATUS MATRIX BANDS SPECIAL CLASS V

For class V restorations (Acc. to Dr. Nyström)
 Stainless steel, Thin 0.040 mm, 12 pcs.

MDD-N1/DS







MDD-N3/DS

MDD-N4A/DS

MDD-N4B/DS

MDD-N5A/DS

MDD-N5B/DS





DENTATUS MO-DO MATRIX BAND

For tight interproximal situations (Acc to Dr. Dyfvermark)
 Stainless Steel, Thin 0.040 mm, 50 pcs.

MDD-HC1 for premolars

MDD-HC2 for molars



DENTATUS MATRIX RETAINERS

MEC

Universal, Straight

MEC-N1

Universal, Angled - Right

MEC-N2


Universal, Angled - Left

MEC-N

Universal, Angled - 2 pcs, Right and Left



MEC



MEC-N

