







Nyström Class V Matrix Band Selection Guide

- Matrix bands for class V cavities near or below the gingival margin.
- Deflects the gingiva to allow access to the cavity
- Provides a dry field of operation

ARTICLE NUMBER	Description	Comments
MDD-N1/DS	“...suitable when cavity is visible from occlusal view”	 N1
MDD-N3/DS	“When the cavity is barely visible”	 N3
MDD-N4A/DS	“When the cavity is quite visible from occlusal view”	 N4A
MDD-N4B/DS	“When the cavity is quite visible from occlusal view” -wider area	 N4B
MDD-N5A/DS	“Especially used on incisors, cuspids or bicuspid for composite in those situations where it is too difficult to apply a rubber dam”	 N5A
MDD-N5B/DS	“Especially used on incisors, cuspids or bicuspid for composite in those situations where it is too difficult to apply a rubber dam” -wider area	 N5B