

JESPER HATT, DDS

In 2003, Dr. Hatt graduated from Aarhus University in Aarhus, Denmark. From 2007 to 2019, he and his wife, Helle Hatt, DDS, had their own fee-for-service practice focusing on comprehensive dentistry. In 2019 they founded AlignerService.com where their team of dentists work as in-house aligner experts online and serve customers in more than 12 different countries. Their services cover all aspects related to clear aligners and, in doing so, they were able to reduce clear aligner therapy chair time by up to 95%.

FOR MORE INFORMATION:

800.323.3136
www.dentatus.com

Complete IPR Solution

An IPR solution with unique features for efficient and safe treatment, as well as increased patient comfort



Find out why Jesper Hatt, DDS, believes that Dentatus' Complete IPR Solution:

- Is more efficient and user-friendly than any other IPR system
- Saves chair time and reduces clinician stress


About 13 years ago, we started using interproximal reduction (IPR) as part of our clear aligner treatment. At the time, IPR was extremely time consuming. We were using metal strips in combination with discs, so the procedure took ages and reduced the profitability of clear aligner therapy dramatically. We needed something more efficient.

In search of something better, I found Dentatus Lamineer Tips. They were ideal for shaping fillings but too small to be used effectively for IPR.

Like any other dentist, I was tempted to try other solutions as they came to the market. For example, several machine-driven reciprocating systems work with short strips. However, in my opinion, all of these systems can be difficult to use when you want to control shaping of the tooth morphology. Especially when compared to my most recent discovery: Dentatus IPR tips and the Profin PDX handpiece.

IPR tips offer several advantages for clinicians, including the following:

- The height of the diamond covered tips is perfect for IPR procedures.
- The reciprocating movement makes it easy to handle stripping and shaping the teeth, and the risk of cutting the patient's mouth is eliminated almost completely.
- The safe and easy handling combined with the intelligent color coding makes it easy to use so IPR can be delegated.



- Whether reducing in the anterior or posterior, Dentatus IPR tips are bliss to use compared to the alternatives.

To my knowledge, Profin PDX is the only handpiece that can switch freely between a free-moving position and a fixed position. In the free-moving position, the IPR tip automatically follows the shape of the teeth and slides into the approximal space. This is helpful when delegating IPR because it minimizes the risk of unintended damage to the teeth. In the fixed position, the tip won't move at all, so it's easy to shape the tooth to the right morphology.

From a cost-benefit perspective, the Dentatus IPR tips and the Profin PDX handpiece are a great add-on for the practice. If you do a lot of clear aligner treatments like we do, the tips will make your IPR procedures easier, faster, and more predictable. They also reduce chair time significantly, which is a benefit for you and your patients.

With the Dentatus IPR solution you can delegate IPR safely to the team with proper training. The result is greater efficiency, increased profitability, and improved patient care.