

DENTATUS POST SYSTEMS - OVERVIEW

Classic Surtex® Posts

The Dentatus Classic Surtex post has a texturized finish, creating two-fold increase of retention compared to machined surfaces. This significantly reduces micro leakage and reduces glare for better aesthetics.


Corresponding Reamers

Helix Classic: RSA #1-6, RSB #1-6, RSC #1-5
Standard Classic: RUA #1-6, RUB #1-6

TITANIUM, SIZE #1-6

GOLD PLATED, SIZE #1-6

STAINLESS STEEL, SIZE #1-6



Luscent Anchors®

Luscent Anchors non-metallic posts transmit light for curing composites within the deep confines of canals providing anchorage and flawless aesthetics. Resin reinforced fiberglass Luscent Anchors fulfill dentistry's preference for very strong, non-metallic posts that are consistent with conventional adhesive techniques and materials to support core foundations and aesthetic restorations. The radioluscent anchors transmit tooth color without being visible under the most adverse light conditions. They also provide a sound core-crown gingival interface radiating natural color through soft and hard tissues surrounding teeth.

Corresponding Reamers

Helix Conical: RSG #3-6, RSH #3-6, RSI #3-5
Standard Conical: RUG #1-6, RUH #4-6

XXXS: LENGTH 15 MM & Ø 1,00 MM


XXS: LENGTH 15 MM & Ø 1,15 MM

XS: LENGTH 15 MM & Ø 1,30 MM

S: LENGTH 15 MM & Ø 1,45 MM

M: LENGTH 17 MM & Ø 1,60 MM

L: LENGTH 19 MM & Ø 1,75 MM



LUSCENT-R Anchors®

LUSCENT-R Anchors offer all the advantages of Luscent Anchors, plus the additional advantage of increased radiopacity levels for clear & crisp imaging. High-quality unidirectional fiber rods are densely distributed and polymerized in a strong resin matrix for superior mechanical strength. The adhesion of the epoxy-resin to the silanized fibers in LUSCENT-R is a key point for excellent fatigue properties.

Corresponding Reamers

Standard Conical: RUG #1-6, RUH #4-6

XXXS: LENGTH 15 MM & Ø 1,00 MM


XXS: LENGTH 15 MM & Ø 1,15 MM

XS: LENGTH 15 MM & Ø 1,30 MM

S: LENGTH 15 MM & Ø 1,45 MM

M: LENGTH 17 MM & Ø 1,60 MM

L: LENGTH 19 MM & Ø 1,75 MM



Twin Luscent Anchors®

Twin Luscent Anchors offer additional advantages to those of Luscent Anchors – improved mechanical retention, both against pulling and rotational forces. It gives a possibility to choose between tapered or parallel post shape, depending on which end is inserted in apical direction.

Corresponding Reamers

Helix Classic: RSA #3-6, RSB #3-6, RSC #3-5
Standard Classic: RUA #3-6, RUB #3-6

XS: LENGTH 16,7 MM & Ø 1,26 MM

S: LENGTH 16,7 MM & Ø 1,40 MM

M: LENGTH 16,7 MM & Ø 1,54 MM





L: LENGTH 16,7 MM & Ø 1,68 MM



CORRESPONDING REAMERS

READ MORE ON PAGES 35-37

POSTS & CORRESPONDING REAMERS

Post system (post sizes)		Corresponding Reamers						
<i>Surtex</i>	<i>Twin Luscent</i>		Helix Classic Reamers				Standard Classic Reamers	
			Short 28 mm	Long 33 mm	X Long 37 mm		Short 28 mm	Long 33 mm
#1	-		RSA-1/3	RSB-1/3	RSC-1/2		RUA-1/3	RUB-1/3
#2	-		RSA-2/3	RSB-2/3	RSC-2/2		RUA-2/3	RUB-2/3
#3	White (XS)		RSA-3/3	RSB-3/3	RSC-3/2		RUA-3/3	RUB-3/3
#4	Yellow (S)		RSA-4/3	RSB-4/3	RSC-4/2		RUA-4/3	RUB-4/3
#5	Red (M)		RSA-5/3	RSB-5/3	RSC-5/2		RUA-5/3	RUB-5/3
#6	Blue (L)	RSA-6/3	RSB-6/3	-	RUA-6/3	RUB-6/3		
<i>Luscent & LUSCENT-R</i>			Conical Reamers				Conical Reamers	
			-	-	-		RUG-1/2	-
	Purple (XXXS)		-	-	-		RUG-2/2	-
	Green (XXS)		-	-	-		RUG-3/2	-
	White (XS)		RSG-3/3	RSH-3/3	RSI-3/2		RUG-4/2	RUH-4/2
	Yellow (S)		RSG-4/3	RSH-4/3	RSI-4/2		RUG-5/2	RUH-5/2
	Red (M)	RSG-5/3	RSH-5/3	RSI-5/2	RUG-6/2	RUH-6/2		
	Blue (L)	RSG-6/3	RSH-6/3	-				

DENTATUS CLASSIC SURTEX® POSTS

PREMIUM DENTAL POSTS WITH PROVEN RELIABILITY

Superior Retention Through Design

The Advantages

- Surtex® surface treatment
- Superior retention
- One-stop dispenser for ease-of-use
- Manufactured in Sweden
- 100% quality control

The original Dentatus Classic Surtex® Posts have become a global standard in dental anchors, and offer dentists a reliable and simple solution for post and core build-up. Our high-quality, pre-fabricated posts are dependable and economic for most endodontic and prosthodontic treatment needs, available in a range of materials including pure titanium, stainless steel and gold plated.


Whatever your crown retention requirements, Dentatus Classic Surtex® Posts has the solution. Available in titanium, stainless steel and gold plated, together with our dedicated reamers they form a state-of-the-art system for post retention in endodontically treated root canals.

Dentatus Classic Surtex® Posts – a continuously developed global standard since 1932.

Square crosscut head for core material adaption

Parallell sides and surface treated for extra retention

Anatomical geometry for tapered roots



CLASSIC SURTEX® POST SIZES

	SHORT (7,8 MM)	MEDIUM (9,3 MM)	LONG (11,8 MM)	EXTRA LONG (14,2 MM)	SUPER LONG (17 MM)
1	Ø 1,05 mm	Ø 1,05 mm	Ø 1,05 mm		
2	Ø 1,20 mm	Ø 1,20 mm	Ø 1,20 mm		
3	Ø 1,35 mm	Ø 1,35 mm	Ø 1,35 mm	Ø 1,35 mm	Ø 1,35 mm
4	Ø 1,50 mm	Ø 1,50 mm	Ø 1,50 mm	Ø 1,50 mm	Ø 1,50 mm
5	Ø 1,65 mm	Ø 1,65 mm	Ø 1,65 mm	Ø 1,65 mm	Ø 1,65 mm
6	Ø 1,80 mm	Ø 1,80 mm	Ø 1,80 mm	Ø 1,80 mm	Ø 1,80 mm

STARTER KITS

CLASSIC SURTEX TITANIUM

Economy Kit

- TST-18** 18 Assorted Posts: 1x S1, S2; 1x M1, M4, M5, M6; 2x M2, M3; 1x L1, L2, L3, L4, L5, L6, 1x EL4, EL6
 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
 1 Hollow Key
 1 Cross Key
 20 Stoppers
 1 One-stop Dispenser



Assortment Kit

- TST-60** 60 Assorted Posts (20 sizes, 3 posts of each)
 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
 1 Hollow Key
 1 Cross Key
 20 Stoppers
 1 One-stop Dispenser




Complete Starter Kit

- TST-B** 120 Assorted Posts (20 sizes, 6 posts of each)
 6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
 1 Hollow Key
 1 Cross Key
 20 Stoppers
 1 One-stop Dispenser



Starter Refill Kit - No Reamers

- TST** 120 Assorted Posts (20 sizes, 6 posts of each)
 1 Hollow Key
 1 Cross Key
 20 Stoppers
 1 One-stop Dispenser



STARTER KITS

CLASSIC SURTEX GOLD PLATED

Economy Kit

- RST-18** 18 Assorted Posts: 1x S1, S2; 1x M1, M4, M5, M6; 2x M2, M3; 1x L1, L2, L3, L4, L5, L6;
1x EL4, EL6
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers
1 One-stop Dispenser



Assortment Kit

- RST-60** 60 Assorted Posts (20 sizes, 3 posts of each)
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers
1 One-stop Dispenser



Complete Starter Kit

- RST-B** 120 Assorted Posts (20 sizes, 6 posts of each)
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers
1 One-stop Dispenser



Starter Refill Kit - No Reamers

- RST** 120 Assorted Posts (20 sizes, 6 posts of each)
1 Hollow Key
1 Cross Key
20 Stoppers
1 One-stop Dispenser



CLASSIC SURTEX STAINLESS STEEL


Economy Kit

- SST-18** 18 Assorted Posts: 1x S1, S2, S3, S4; 1x M1, M2, M3, M4; 1x L3, L4, L5, L6; 1x EL3, EL4, EL5;
1x SL3, SL4, SL5
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers
1 One-stop Dispenser




Assortment Kit

- SST-60** 60 Assorted Posts (20 sizes, 3 posts of each)
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers
1 One-stop Dispenser



Complete Starter Kit

- SST-B** 120 Assorted Posts (20 sizes, 6 posts of each)
6 Assorted Standard Classic Reamers (1 RUA of each size, length 28 mm)
1 Hollow Key
1 Cross Key
20 Stoppers
1 One-stop Dispenser



Starter Refill Kit - No Reamers

- SST** 120 Assorted Posts (20 sizes, 6 posts of each)
1 Hollow Key
1 Cross Key
20 Stoppers
1 One-stop Dispenser



REFILLS

CLASSIC SURTEX TITANIUM

Short, 7.8 mm, 15 pcs/pkg

TST-S1	Diameter 1.05 mm
TST-S2	Diameter 1.20 mm
TST-S3	Diameter 1.35 mm
TST-S4	Diameter 1.50 mm
TST-S5	Diameter 1.65 mm
TST-S6	Diameter 1.80 mm



Medium, 9.3 mm, 15 pcs/pkg

TST-M1	Diameter 1.05 mm
TST-M2	Diameter 1.20 mm
TST-M3	Diameter 1.35 mm
TST-M4	Diameter 1.50 mm
TST-M5	Diameter 1.65 mm
TST-M6	Diameter 1.80 mm



Long, 11.8 mm, 15 pcs/pkg

TST-L1	Diameter 1.05 mm
TST-L2	Diameter 1.20 mm
TST-L3	Diameter 1.35 mm
TST-L4	Diameter 1.50 mm
TST-L5	Diameter 1.65 mm
TST-L6	Diameter 1.80 mm



Extra Long, 14.2 mm, 15 pcs/pkg

TST-EL3	Diameter 1.35 mm
TST-EL4	Diameter 1.50 mm
TST-EL5	Diameter 1.65 mm
TST-EL6	Diameter 1.80 mm



Super Long, 17 mm, 15 pcs/pkg

TST-SL3	Diameter 1.35 mm
TST-SL4	Diameter 1.50 mm
TST-SL5	Diameter 1.65 mm
TST-SL6	Diameter 1.80 mm



REFILLS

CLASSIC SURTEX GOLD PLATED

Short, 7.8 mm, 15 pcs/pkg

RST-S1	Diameter 1.05 mm
RST-S2	Diameter 1.20 mm
RST-S3	Diameter 1.35 mm
RST-S4	Diameter 1.50 mm
RST-S5	Diameter 1.65 mm
RST-S6	Diameter 1.80 mm



Medium, 9.3 mm, 15 pcs/pkg

RST-M1	Diameter 1.05 mm
RST-M2	Diameter 1.20 mm
RST-M3	Diameter 1.35 mm
RST-M4	Diameter 1.50 mm
RST-M5	Diameter 1.65 mm
RST-M6	Diameter 1.80 mm



Long, 11.8 mm, 15 pcs/pkg

RST-L1	Diameter 1.05 mm
RST-L2	Diameter 1.20 mm
RST-L3	Diameter 1.35 mm
RST-L4	Diameter 1.50 mm
RST-L5	Diameter 1.65 mm
RST-L6	Diameter 1.80 mm



Extra Long, 14.2 mm, 15 pcs/pkg

RST-EL3	Diameter 1.35 mm
RST-EL4	Diameter 1.50 mm
RST-EL5	Diameter 1.65 mm
RST-EL6	Diameter 1.80 mm



Super Long, 17 mm, 15 pcs/pkg

RST-SL3	Diameter 1.35 mm
RST-SL4	Diameter 1.50 mm
RST-SL5	Diameter 1.65 mm
RST-SL6	Diameter 1.80 mm



REFILLS

CLASSIC SURTEX STAINLESS STEEL

Short, 7.8 mm, 12 pcs/pkg

SST-S1	Diameter 1.05 mm
SST-S2	Diameter 1.20 mm
SST-S3	Diameter 1.35 mm
SST-S4	Diameter 1.50 mm
SST-S5	Diameter 1.65 mm
SST-S6	Diameter 1.80 mm



Medium, 9.3 mm, 12 pcs/pkg

SST-M1	Diameter 1.05 mm
SST-M2	Diameter 1.20 mm
SST-M3	Diameter 1.35 mm
SST-M4	Diameter 1.50 mm
SST-M5	Diameter 1.65 mm
SST-M6	Diameter 1.80 mm



Long, 11.8 mm, 12 pcs/pkg

SST-L1	Diameter 1.05 mm
SST-L2	Diameter 1.20 mm
SST-L3	Diameter 1.35 mm
SST-L4	Diameter 1.50 mm
SST-L5	Diameter 1.65 mm
SST-L6	Diameter 1.80 mm



Extra Long, 14.2 mm, 12 pcs/pkg

SST-EL3	Diameter 1.35 mm
SST-EL4	Diameter 1.50 mm
SST-EL5	Diameter 1.65 mm
SST-EL6	Diameter 1.80 mm



Super Long, 17 mm, 12 pcs/pkg

SST-SL3	Diameter 1.35 mm
SST-SL4	Diameter 1.50 mm
SST-SL5	Diameter 1.65 mm
SST-SL6	Diameter 1.80 mm



ACCESSORIES

FOR DENTATUS METAL POSTS

NLH	1 Hollow Key
NLK	1 Cross Key

