DENTATUS IPR TIPS FOR PROFIN

FOR ALIGN ORTHODONTICS

Your Complete IPR Solution

The Advantages:

- Efficient and safe treatment
- Increased patient comfort
- Gentle to the gingiva
- One-hand operation

The solution consists of Interproximal Reduction (IPR) tips together with the reciprocating Profin PDX handpiece.



Profin PDX with IPR Tip




IPR in Align Orthodontics



Clear Aligner Treatment

Calibrated blade thickness from 0.1 mm to 0.5 mm Long lasting tips Autoclavable



Art.nr. LTA-P/SK2



Control and Performance Excellence with Profin's unique features

The Profin PDX handpiece has **unique features** such as reciprocating with a 1.2 mm axial stroke and self-steering tips. Alternatively the head has the possibility to lock the reciprocating tips in 8 different positions. Profin IPR's controlled reduction avoids rounding teeth thus preserving the natural anatomy of the tooth.IPR tips are color-coded to ensure time and performance excellence.

Product Description

You may choose to buy an IPR Starter Kit (**Art.nr. LTA-P/SK2**), including 5 tips and a stand or your chosen tips separately in a IPR refill pack of 2 pcs.

IPR Starter Kit



LTA-P/SK2

LTA-P4: Green	Double sided. Blade thickness 0.4 mm. Grit size 50µ.
LTA-P5 : Red	Double sided. Blade thickness 0.5 mm. Grit size 50µ.

Accessory Tips

LTA-P25R: Violet	Safe-sided – abrasive on the right side. Blade thickness 0.25 mm. Grit size 50µ.
LTA-P25L: Violet	Safe-sided – abrasive on the left side. Blade thickness 0.25 mm. Grit size 50µ.



LTA-P3/2

IPR Refills Kits

All IPR Tips are color coded and have the grit size equivalent to the description below. Their colorcoding is strictly related to blade thickness.

LTA-P1: Black	Double sided. Blade thickness 0.1 mm. Grit size 15µ.
LTA-P2: White	Double sided. Blade thickness 0.2 mm. Grit size 40µ.
LTA-P3: Yellow	Double sided. Blade thickness 0.3 mm. Grit size 50µ.

Safety Excellence

As safety is the ultimate concern for the dentist, we have developed two variants of our IPR tips (see product description below for details). IPR tips are either coated double-sided with diamond or safe-sided to ensure safety to adjacent teeth.



Dentatus AB www.dentatus.com SE: +46 8 546 509 00 info@dentatus.se US: +1 212 481 1010 dentatus@dentatus.com

