

**Surtex**<sup>®</sup>  
STERILE  
endodontic posts



# Bacteria not included.

**There's never been a better time to switch to Dentatus Posts.  
Announcing Surtex Sterile, the first and only pre-sterile post system—  
providing safety and confidence to clinician, staff and patients alike.**

## **PREMIUM ENDODONTIC POSTS WITH PROVEN RELIABILITY Superior Retention Through Design**

### **The post**

- Surtex<sup>®</sup> surface treatment = superior retention
- Square crosscut head for core material adaption
- Parallel sided and surface treated for extra retention
- Anatomical geometry for tapered roots
- Threaded for easy retrieval
- Available in Gold-Plated or Titanium

### **The quality**

- Manufactured in Sweden
- 100% quality control

### **The sterile packaging**

- Individual cavities for each sterile post
- Perforated & easy to separate
- Color coded and easy to organize
- Sterile posts delivery for maximum safety



Available in Titanium and Gold Plated.

**(800) 323-3136**  
**www.dentatus.com**



# Surtex® STERILE endodontic posts

**Dentatus** is a global expert in the endodontic post market, offering dentists consistently reliable solutions for post and core build-ups. Dentatus revolutionized dentistry in 1932 manufacturing the first prefabricated endodontic post and they're doing it again with the first commercially available pre-sterile post.



GOLD-PLATED

**Surtex® Sterile Posts** are pre-abraded and pre-sterilized saving you time and ensuring compliance in infection control protocols. Packaged in individual blister packs of 5 posts in Titanium and Gold Plated, Surtex Sterile are easy to organize and color coded for quick reference. Learn more at [Dentatus](http://Dentatus.com) | [Surtex Sterile](http://SurtexSterile.com).

## Ordering Information

LENGTH	DIAMETER	GOLD PLATED	TITANIUM
Short 7.8mm	1.05 mm	RSTS-S1	TSTS-S1
	1.05 mm	RSTS-S2	TSTS-S2
	1.35 mm	RSTS-S3	TSTS-S3
	1.50 mm	RSTS-S4	TSTS-S4
	1.65 mm	RSTS-S5	TSTS-S5
	1.80 mm	RSTS-S6	TSTS-S6
Medium 9.3mm	1.05 mm	RSTS-M1	TSTS-M1
	1.05 mm	RSTS-M2	TSTS-M2
	1.35 mm	RSTS-M3	TSTS-M3
	1.50 mm	RSTS-M4	TSTS-M4
	1.65 mm	RSTS-M5	TSTS-M5
	1.80 mm	RSTS-M6	TSTS-M6
Long 11.8mm	1.05 mm	RSTS-L1	TSTS-L1
	1.05 mm	RSTS-L2	TSTS-L2
	1.35 mm	RSTS-L3	TSTS-L3
	1.50 mm	RSTS-L4	TSTS-L4
	1.65 mm	RSTS-L5	TSTS-L5
	1.80 mm	RSTS-L6	TSTS-L6
Extra Long 14.2mm	1.35 mm	RSTS-EL3	TSTS-EL3
	1.50 mm	RSTS-EL4	TSTS-EL4
	1.65 mm	RSTS-EL5	TSTS-EL5
	1.80 mm	RSTS-EL6	TSTS-EL6
Super Long 17mm	1.35 mm	RSTS-SL3	TSTS-SL3
	1.50 mm	RSTS-SL4	TSTS-SL4
	1.65 mm	RSTS-SL5	TSTS-SL5
	1.80 mm	RSTS-SL6	TSTS-SL6

## SURTEX STERILE ENDODONTIC POSTS 5/PK

### NEW TO SURTEX?

Be sure to include the instrumentation kit in your starting order [RUA].

**RETURNING CUSTOMERS:** Thank you, you're all set! Surtex Sterile works with Surtex instruments, just select your post sizes.

**We recommend stocking the 5 most common sizes: M2, M3, L2, L3 & L4**

### HOW TO ORDER THE CORRECT POST:

#### SELECT THE MATERIAL

RSTS= Gold-Plated | TSTS= Titanium

#### SELECT THE LENGTH

Length (L)

S Short 7.8mm

M Medium 9.3mm

L Long 11.8mm

EL Extra Long 14.2mm

SL Super Long 17mm

#### SELECT THE DIAMETER

Diameter (Ø)

1 1.05 mm ●

2 1.20 mm ●

3 1.35 mm ●

4 1.50 mm ●

5 1.65 mm ●

6 1.80 mm ●

**BE SURE TO INCLUDE THE INSTRUMENTATION KIT IN YOUR STARTING ORDER**

**Instrumentation Kit [RUA]**  
Contains: 6 Dentatus Classic Reamers [RUA] 28mm long and size matched to the 6 Surtex Sterile post diameters.

(800) 323-3136  
[www.dentatus.com](http://www.dentatus.com)

