

# DENTATUS MATRIX BANDS

THE ORIGINAL & SPECIAL CLASS V

## Highly Regarded by Users All Over the World

### The Advantages


- Swedish stainless steel
- Recognized and used globally
- Original Nyström retainers – no copies

Dentatus' Matrix Bands are manufactured from the highest quality stainless steel. The standard shape bands are available in two thicknesses: 0.050 and 0.040 mm.

The 0.040 bands are made of high quality stainless steel, and the thicker bands, 0.050, are made from stainless steel with special dead soft quality. The choice of steel quality gives the bands high strength and good adjustability. The thinnest bands provide very good prerequisites for well contoured and defined interproximal contacts.

Dentatus quality has proven to be highly regarded by users all over the world for a long period of time.

All Dentatus Matrix Bands have perfectly smooth edges.




# Original Dentatus Matrix Bands

Stainless Steel, Thin 0.040 mm, 12 pcs/pkg

Stainless Dead Soft Steel, 0.050 mm,  
12 pcs/pkg and 100 pcs/pkg

MDD-1RF MDD-2RF MDD-3RF  
MDD-4RF MDD-5RF MDD-6RF


MDD-1 MDD-2 MDD-3  
MDD-4 MDD-5 MDD-6



# Dentatus Matrix Bands Special Class V

For Class V Restorations, Stainless Steel, Thin 0.040 mm (acc to Dr Nyström), 12 pcs/pkg

MDD-N1/DS MDD-N3/DS MDD-N4A/DS MDD-N4B/DS MDD-N5A/DS MDD-N5B/DS



Dentatus AB  
www.dentatus.com  
SE: +46 8 546 509 00  
info@dentatus.se  
US: +1 212 481 1010  
dentatus@dentatus.com

